

Timing Filter Amplifier



- Timing with germanium detectors
- Energy spectroscopy at ultra-high count rates
- Selectable filter for pulse shaping
- Signal-to-noise ratio optimization
- Continuously adjustable gain X2 to X250
- Pole-zero cancellation
- DC-coupling

The Model 474 Timing Filter Amplifier is especially designed to shape pulses and permit optimizing the signal-to-noise ratio for timing measurements. The Model 474 is particularly suited for use with an ORTEC Constant-Fraction Timing Discriminator in timing applications with germanium or silicon charged-particle detectors ([Fig. 1](#)). It derives its input signal directly from the preamplifier output. The timing spectrum in [Fig. 2](#) illustrates the performance obtainable with the Model 474 shaping the germanium detector pulses before they are furnished to the discriminator. Tables [1](#) and [2](#) give typical performance data for various ORTEC germanium detectors.

The fast rise time, high output drive, and high gain capabilities of the Model 474 make it useful for other applications, such as timing with systems utilizing low-gain photomultiplier tubes. In addition, the pole-zero cancellation network, the dc-coupling, and the time-invariant baseline restorer permit energy spectroscopy with scintillation detectors and Si charged-particle detectors at ultra-high count rates.

A wide variety of input pulse shapes can be filtered as required for optimum signal processing. The Model 474 combines continuously adjustable gain (X2 to X250) with separately selectable Integrate (τ_i) and Differentiate (τ_d) time constants for proper pulse shaping, making this unit an important asset for time measurement.

PERFORMANCE

Input Amplitude Range 0 to ± 1 V signal, 0 to ± 5 V dc offset; maximum input ± 5 V total.

Output Amplitude Range 0 to ± 5 V with a 50- Ω load.

Noise For maximum gain, rms noise referred to the input is ≤ 10 μ V with $\tau_i = \tau_d = 200$ ns or ≤ 50 μ V with filter out; measured using a Hewlett-Packard 3400A true rms meter.

Rise Time ≤ 10 ns with filter out or $\sim 2.2 \tau_i$ for other selections.

Nonlinearity $\leq \pm 0.05\%$ at midband frequency over ± 5 V range.

Temperature Instability

DC Level $\leq \pm 25 \mu$ V/ $^{\circ}$ C referred to the output.

Gain $\leq \pm 0.06\%$ / $^{\circ}$ C.

Specifications over 0 to 50 $^{\circ}$ C range.

CONTROLS

COARSE GAIN Front-panel 6-position switch for selection of X1, X2, X4, X6, X10, or X20.

FINE GAIN Front-panel single-turn potentiometer, continuous from X2 to X12.5.

PZ ADJ Front-panel screwdriver adjustment to compensate for the preamplifier decay time constant.

Time Constant Two 6-position switched on front panel:

INTEGRATE RC time constants: OUT (equivalent to 4 ns), 20, 50, 100, 200, and 500 ns.

DIFF RC time constants: OUT (equivalent to 0.2 ms), 20, 50, 100, 200, and 500 ns.

NOTE: With Differentiate and Integrate in the OUT position, the passband is 1 kHz to 35 MHz.

N30526A	73%	2.03	34	5.4	19.4	3.7	8.8
P20171	81%	1.97	34	5.5	27.0	4.7	13.8
N20366A	88%	2.34	36	5.8	21.2	5.5	16.4

Fig. 1 Simplified Timing System.




Fig. 2 Timing Spectrum for a Narrow Synamic Range Using a Germanium Detector. (Resolution values are given in Table 1.)

22Na
 Start: KL236 (1x1), RCA 8575
 Photomultiplier Tube
 Stop: Ge Coax, 12.5%, 62.3 cc
 1.1:1 Dynamic Range
 473A Discriminator

