

PRELIMINARY

INSTRUCTION MANUAL

TC 178/TC 178P

QUAD PREAMP/QUAD PULSER

NSCL-ELECTRONIC

WARRANTY

TENNELEC, INC. warrants that the products or components manufactured by it shall be free from defects in material or workmanship for a period of two years from **the date of delivery** to purchaser. If such product or component is determined to be defective by TENNELEC, its sole warranty obligation shall be limited to either replacing or repairing such defective product or component **or allowing credit** therefor, at **TENNELEC's** option. Such warranty is further conditioned upon the purchaser's giving prompt notice of any such defect and satisfactory proof thereof to **TENNELEC's** customer service manager, thereafter upon TENNELEC's approval, the purchaser shall return such defective product or component to **TENNELEC's** factory at Oak Ridge, Tennessee, all transportation charges prepaid. TENNELEC shall be responsible only for transportation charges incurred in returning **such** product or component to purchaser. All customs, brokerage and duty charges shall be at the expense of the purchaser. Damage in transit due to inadequate packaging will be repaired at **purchaser's** expense. Any repairs **or** replacements by the purchaser without **TENNELEC's** approval, any willful abuse or any evidence that the product or component was not properly used and maintained, would automatically void this warranty.

TENNELEC makes no warranty whatsoever in respect to products or components not manufactured by it but instead the applicable warranties, if any, of the respective manufacturers thereof shall **apply**. Likewise fuses, batteries and input transistors in ultra low-noise amplifiers are specifically excluded from this warranty.

THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, STATUTORY OR OTHERWISE, INCLUDING WARRANTY OF MERCHANTABILITY AND FITNESS.

TABLE OF CONTENTS

Section	Page
1.0 INTRODUCTION.	1
1.1 TC 178 QUAD PREAMPLIFIER	1
1.2 TC 178P QUAD PULSER	1
2.0 SPECIFICATIONS	2
2.1 TC178	2
2.1.1 PERFORMANCE.	2
2.1.2 CONTROLS	6
2.1.2A FRONT-PANEL CONTROLS	6
2.1.2B INTERNAL CONTROLS	6
2.1.3 CONNECTORS	7
2.1.3A FRONT PANEL	7
2.1.3B REAR PANEL	7
2.1.4 INDICATORS	7
2.1.5 POWER REQUIREMENTS	8
2.1.6 OTHER INFORMATION.	8
2.2 TC 178P	8
2.2.1 PERFORMANCE	8
2.2.2 CONTROLS	9
2.2.3 CONNECTORS
2.2.3A FRONT PANEL	9
2.2.3B REAR PANEL	10
2.2.4 POWER REQUIREMENTS	10
2.2.5 OTHER INFORMATION	10

Section	Page
3.0 INSTALLATION	11
3.1 TC178	11
3.1.1 POWER	1
3.1.2 CONNECTIONS	11
3.1.2A DETECTOR CONNECTION	11
3.1.2B BIAS SUPPLY CONNECTION	12
3.1.2C ENERGY OUTPUT CONNECTION	12
3.1.2D TEST IN CONNECTION	13
3.1.3 GENERAL PRECAUTIONS	14
3.1.3A FET PROTECTION	14
3.1.3B APPLYING BIAS VOLTAGE	14
3.1.3C DETECTOR BIAS NETWORK	15
3.1.3D DETECTOR BIAS LEVEL	15
3.2 TC 178P	16
3.2.1 POWER	16
3.2.2 CONNECTIONS	17
3.2.2A PULSER OUT	17
3.2.2B SYNC OUT	18
3.2.3C EXT IN	18
4.0 OPERATING PROCEDURES	18
4.1 FIRST-TIME OPERATION	18
4.1.1 EQUIPMENT REQUIRED	18
4.1.2 INITIAL TEST SETUP	19
TABLE OF CONTENT;	19

TABLE OF CONTENTS

Section	Page
4.1.3 PULSER CALIBRATION	22
4.2 ADJUSTMENTS	23
4.2.1 TC178	23
4.2.1A DC OFFSET VOLTAGE	23
4.2.1B RISETIME	23
4.2.1C PULSER CALIBRATION	25
4.2.2 TC 178P	26
5.0 NOISE	26
5.2 NOISE MEASUREMENTS	30
6.0 COUNT RATE EFFECTS	31
6.1 RESOLUTION	31
6.2 NONLINEARITY	34
7.0 PREAMPLIFIER MODIFICATIONS	34
7.1 CHANGING PREAMPLIFIER SENSITIVITY	34
7.2 FRONT-END MODIFICATIONS	35
7.3 MODIFICATION OF FET DRAIN CURRENT	35
7.4 ELECTRONIC NOISE REDUCTION	35
8.0 SHIPPING DAMAGE	36
9.0 SERVICING	36
10.0 WARRANTY	36

TABLE OF FIGURES

Figure	Page
2.1 Electronic Noise vs Equivalent Detector Capacitance (Peaking Time of 4.0 usec)	3
2.2 Electronic Noise vs Equivalent Detector Capacitance (Peaking Time of 1.0 usec)	3
2.3 Typical Risetime vs Equivalent Detector Capacitance	4
2.4 Typical Transient Response A) 0.1 GeV , B) 1.0 GeV , C) 10 GeV	5
2.5 Typical Sync Out Signal	9
2.6 Typical Pulser Out Signal	10
3.1 External Test Input Terminator	13
4.1 Initial Test Setup	20
4.2 Typical TC 178 Output, 0.1 GeV Range	20
4.3 Typical TC 178 Output, 1.0 GeV Range	21
4.4 Typical TC 178 Output, 10.0 GeV Range,	21
5.1 Electronic Noise vs Equivalent-Detector Capacitance	27
5.2 Electronic Noise vs Equivalent Detector Capacitance	28
5.3 Electronic Noise vs Equivalent Detector Capacitance	29
6.1 The Shape of a Typical Spectral Line (Gaussian).	31
6.2 Pulse Shape at the Output of a Preamplifier	32
6.3 The Waveform at the Output of a Preamplifier with Pulser Applied to the Input in Rapid Succession	33
TC 178 SCHEMATIC DIAGRAM (Sheet 1 of 2	38
TC 178 SCHEMATIC DIAGRAM (Sheet 2 of 2)	39
TC 178P SCHEMATIC DIAGRAM (Sheet 1 of 1)	40

1.0

The TENNELEC TC 178 and TC 178P are designed to operate together as a pair. The TC 178 is a four channel, ac-coupled charge sensitive preamplifier designed for use with detectors having an equivalent capacitance of from 0pF to greater than 1000pF. The TC 178P is a four channel pulser used to drive the TC 178 internal pulser and provide amplitude and rate control.

1.1 TC 178 QUAD PREAMPLIFIER

The TC 178 Quad Preamplifier contains four independent charge-sensitive preamplifiers. Each channel has three independent energy ranges selected via front-panel mounted toggle switches. The internal pulser polarity is selected independently for each channel via toggle switches accessible through the bottom of the TC 178 case. The pulser calibration for each channel is accessible through the bottom of the preamplifier case. The pulser calibration is independent for each channel and energy range. Visual indication of pulse processing for each channel is provided by LED indicators. A preamplifier power LED indicator illuminates when supply voltage levels are within acceptable limits and extinguishes if any of the supply voltages fall outside these limits.

Each preamplifier channel accepts detector bias up to 1.0kV through an SHV connector. The detector bias network has a total resistance of 3.0 megohm, of which, 1.0 megohm is the detector load resistor. The detector signal is coupled to the preamplifier input via SHV connectors.

The risetime of each preamplifier is independently adjustable for optimizing the preamplifier response with various detectors. The preamplifier dc-off set adjustment is available to the user but will generally only need adjustment if the input FET is replaced.

1.2 TC 178P QUAD PULSER

The TC 178P Quad Pulser is a single-wide NIM that is specifically dedicated for use with the TC 178 Quad Preamplifier. The TC 178P provides a dc reference level and a trigger pulse to the TC 178 TEST IN connector. The dc reference is independent for each channel and is used to determine the level of the pulse.

generated by the internal pulser of the TC 178. The rate of operation is common to all four channels and is continuously variable from 1 Hz to 11 kHz through four ranges. An **external input** is provided to **allow** inputs from random pulsers or rates up to approximately 100 kHz.

The TC 178P, when used with the TC 178, can be calibrated to read directly in terms of equivalent energy. The calibration feature allows for both testing and calibration of nuclear counting system electronics.

2.0 **SPECIFICATIONS**

2.1 TC 178 (per channel except where noted).

2.1.1 **PERFORMANCE***

NOISE (FWHM referred to a **silicon** detector with $W = 3.6 \text{ eV}$ per electronhole pair).

DETECTOR CAPAC. (pF)	NOISE (keV) FWHM (1)								
	ENERGY RANGE								
	0.1 GeV			1.0 GeV			10 GeV		
	$t_p=4\text{usec}$		$t_p=1\text{usec}$	$t_p=4\text{usec}$		$t_p=1\text{usec}$	$t_p=4\text{usec}$		$t_p=1\text{usec}$
	TYP	MAX	TYP	TYP	MAX	TYP	TYP	MAX	TYP
0	11.6	25.0	11.0	14.6	30.0	16.0	71	150	112
10	11.8		11.7	14.8		16.6	71		112
20	12.0		12.0	15.0		16.9	71		113
50	12.6		13.3	15.6		18.0	72		113
100	14.0	35.0	15.7	16.7	40.0	20.0	73	160	114
200	16.8		21.4	19.6		25.0	76		121
500	27.0		38.5	29.2		41.0	85		132
1000	43.7	110.0	68.0	47.0	120.0	70.0	100	250	153

TYP INTERCEPT 11.0keV TYP INTERCEPT 16.0keV TYP INTERCEPT 71.0keV
 TYP SLOPE 30 eV/pF TYP SLOPE 30 eV/pF TYP SLOPE 30 eV/pF

(1) t_p refers to the peaking time in usec. The peaking time is a measure of the time from the 1% level to the 100% level of the shaped pulse. The peaking time is approximately two times the shaping time. The noise measurements were made using a TC 243 shaping amplifier in the triangular mode.

*Where appropriate, all performance specifications given are for a 4 usec peaking (2 usec shaping) time.

See Figure 2.1 and Figure 2.2 for graphs of TC 178 Electronic Noise vs Equivalent Detector Capacitance at peaking times of 4 usec and 1 sec.




Figure 2.1 Electronic Noise vs Equivalent Detector Capacitance (Peaking Time of 4.0 usec)




Figure 2.2 Electronic Noise vs Equivalent Detector Capacitance (Peaking Time of 1.0 usec)

RISETIME (2)

DETECTOR CAPAC. (pF)	RISETIME nsec (10% - 90%)					
	ENERGY RANGE					
	0.1 Gev		1.0 Gev		10 Gev	
	TYP	MAX	TYP	MAX	TYP	MAX
0	17	25	14	20	10	15
10	17		14		10	
20	19		14		10	
50	21		15		11	
100	25		21		13	
200	38	50	48	25	19	20
500	80					
1000	160		100		80	

(2) Based on a full scale equivalent input, **risetime** adjustment optimized at each measurement, E OUT terminated in 50 ohms, measurement made with an external test-input capacitor, and internal protection network installed.

See Figure 2.3 for a graph of the TC 178 **Risetime** as a function of **equivalent detector capacitance**. The typical transient response of the TC 178 is shown in Figure 2.4.





Figure 2.3 Typical **Risetime vs** Equivalent Detector Capacitance



A

B

C

Figure 2.4 Typical Transient Response
A) 0.1 GeV, B) 1.0 GeV, C) 10 GeV

CONVERSION GAIN (Si)

RANGE

0 - 0.1 GeV

20 mV/MeV

0 - 1.0 GeV

2.0 mV/MeV

0 - 10 GeV

0.2 mV/MeV

CONVERSION GAIN TEMPERATURE
INSTABILITY

RANGE

0 - 1.0 GeV

$\leq \pm 100$ ppm/°C

0 - 1.0 GeV

$\leq \pm 75$ ppm/°C

0 - 10 GeV

$\leq \pm 50$ ppm/°C

OPERATING TEMPERATURE RANGE

0 to 50°C

INTEGRAL NONLINEARITY

0 to ± 10 V

1.05% MAX., typically 10.02%

0 to ± 2 V

$\leq 0.02\%$ MAX., typically 10.005%

DIFFERENTIAL NONLINEARITY

0 to ± 10 V

Typically 10.03%

DYNAMIC INPUT
CAPACITANCE

RANGE

0 - 0.1 GeV

0 - 1.0 GeV

0 - 10 GeV

$\geq 10\text{kpF}$; typically $>15\text{kpF}$

$\geq 20\text{kpF}$; typically $>60\text{kpF}$

$\geq 30\text{kpF}$; typically, $>85\text{kpF}$

MAXIMUM ENERGY

RANGE

0 - 0.1 GeV

0 - 1.0 GeV

0 - 10 GeV

$>500\text{ MeV (Si)}$

$>5\text{ GeV (Si)}$

$>50\text{ GeV (Si)}$

COUNT RATE CAPABILITY

(5% of pulses in
non-linear range)

RANGE

0 - 0.1 GeV

0 - 1.0 GeV

0 - 10 GeV

$>1 \times 10^5$ cps @ 100 MeV (Si)

$>1 \times 10^5$ cps @ 1.0 GeV (Si)

$>1 \times 10^5$ cps @ 10 GeV (Si)

DECAY TIME CONSTANT
(NOMINAL)

50 usec

DETECTOR LOAD RESISTOR

1 megohm

DETECTOR BIAS CAPABILITY

$\pm 1,000\text{V}$ Maximum

SIGNAL POLARITY

INWT

E CDT

TEST IN

POSITIVE OR NEGATIVE

INVERSE OF INPUT

POSITIVE (dc level plus
trigger pulse from TC 178P)

2.1.2 CONTROLS

2.1.2A FRONT-PANEL CONTROLS

RANGE Locking toggle switch for selection of energy
range, 0.1, 1.0, or 10 GeV full scale.

2.1.2B INTERNAL CONTROLS

DC-OFFSET Screwdriver-adjustable multiturn
potentiometer (accessible through the top of the
preamplifier case) used to adjust input FET drain
current and dc offset voltage. NO user adjustment
normally required unless input FET is replaced. One
per channel.

RISETIME Screwdriver adjustable single-turn potentiometer (accessible through the top of the preamplifier case) used to optimize the preamplifier transient response for each detector preamplifier **combination**. One per channel.

PULSER POLARITY Two-position toggle switch selects either positive (POS) or **negative (NEG)** pulser polarity. One per channel, accessible through bottom of case.

PULSER CALIBRATION Three variable capacitors, one per **channel, accessible** through bottom of case, used to **calibrate pulser** amplitude for each range.

2.1.3 CONNECTORS

2.1.3A FRONT PANEL

INPUT SBV (AMP 51494-2) connector accepts positive or negative **charge** from the detector and applies high voltage to the detector. One per channel.

2.1.3B REAR PANEL

H.V. IN SHV (AMP 51494-2) connector accepts up to **±1000V** and **applies** it to the detector through the INPUT connector. **One** per channel.

TEST IN BNC (UG-1094/U) connector accepts pulser level and trigger **signals** from TC 178P. One per channel.

E OUT BNC (UG-1094/U) connector provides energy output signal, **dc-coupled, $Z_0 = 50$ ohm $\pm 1\%$** , dc offset approximately -100 **mV**. One per channel.

POWER **9-pin** male, **Amphenol 17-20090** or equivalent. One per TC 178.

2.1.4 INDICATORS

E-OUT LED is illuminated when the E-OUT level exceeds approximately 2% of **full** scale range (2 **MeV**, 20 **MeV** or 200 **MeV**).

POWER ON LED is illuminated when supply voltage levels are within acceptable limits} extinguishes if any **of** the supply voltages fall outside of these limits.

2.1.5 **POWER REQUIREMENTS +24V @ 160 mA +12V @ 210 mA**
-24V @ 130 mA -12V @ 0 mA

2.1.6 OTHER INFORMATION

WEIGHT (SHIPPING) 6.0 lbs (2.7 kg)
(NET) 4.0 lbs (1.8 kg)

DIMENSIONS (L X W X H) 8.5 x 5 x 3 inches
(21.6 X 12.7 X 7.6 cm)

INSTRUCTION MANUAL One provided with each **instrument ordered.**

ACCESSORY One **TENNELEC NC-PAC-10, 10 ft.** preamplifier power cable provided with each **preamplifier ordered; Amphenol 17-20090 to Amphenol 17-10090 connectors.**

2.2 **TC 178P**

2.2.1 PERFORMANCE

REPETITION RATE 1 Hz to 11 **kHz** in 4 ranges; continuously variable for each range

RATE INACCURACY $\leq \pm 2\%$ of setting, typically $< \pm 1\%$

RATE TEMPERATURE INSTABILITY $\leq \pm 200 \text{ ppm}/^{\circ}\text{C}$, typically $\leq \pm 50 \text{ ppm}/^{\circ}\text{C}$ referenced to **+24V**

RATE INSTABILITY $\leq 4.25\%/V$ for variations of the **+24V** NIM supply voltage; independent of **115V** line voltage variations

PULSE AMPLITUDE RANGE 0 to **+10.0V** dc

PULSE AMPLITUDE INACCURACY $< \pm [25 \text{ mV} + (V_o \text{ DIAL} \times 2.8 \text{ mV})]$ typically $\leq \pm 10 \text{ mV}$

PULSE AMPLITUDE TEMPERATURE INSTABILITY $V_o < 0.4V$ $\leq \pm 50 \text{ ppm}/^{\circ}\text{C}$; typically $\leq \pm 20 \text{ ppm}/^{\circ}\text{C}$

$V_o \geq 0.4V$ $\leq \pm 40 \text{ ppm}/^{\circ}\text{C}$; typically $\leq \pm 20 \text{ ppm}/^{\circ}\text{C}$

PULSE **AMPLITUDE** ≤ 0.1 mV/V for variations of
INSTABILITY the +24V NIM supply voltage;
independent of 115V line
variations.

2.2.2 CONTROLS

RATE Front-panel, six-position rotary switch selects pulse repetition rates of 1, 10, 100, **1kHz**, OFF, or **external** (EXT) for all four channels.

MULT Front-panel 10-turn potentiometer provides multiplier for **repetition rate** for all four channels! factors are X1.0 to X11, continuously variable.

PULSE AMPLITUDE Front-panel 10-turn potentiometer provides **adjustment of output dc** level for each channel from 0 to **+10V** open circuit; One per channel.

ON-OFF **Two-position toggle** switch to enable or **disable** the **respective pulser** channel. One per **channel**.

2.2.3 CONNECTORS

2.2.3A FRONT PANEL

EXT IN Front-panel BNC (UG-1094/U) connector, **dc-**coupled, accepts positive NIM input (**+5V, 500 nsec** wide) to provide pulse repetition rates **up to 100 kHz** for all four channels. The EXT IN signal minimum usable width **is** typically **≤ 100 nsec**. The, **EXT IN** connector is duplicated on the Rear Panel.

SYNC OUT Front-panel BNC (UG-1094/U) **connector** provides de-coupled, 0 to **+5V** signal open circuit (0 to **+2.5V** terminated into **50 ohms**). Signal width is typically 60 naec. The SYNC OUT connector is duplicated on the Rear Panel. See Figure 2.5 for a typical SYNC OUT signal.




Figure 2.5 Typical Sync Out Signal

2.2.3B REAR PANEL

PULSER OUT Rear-panel BNC (Amphenol **31-010**, insulated **BNC**), dc-coupled, provides the TC **178** trigger pulse (TTL level) and the TC 178 level control, **0** to **+10.0V** dc open circuit. One per channel. See Figure 2.6 for a typical PULSER OUT signal.




Figure 2.6 Typical Pulser Out Signal

EXT IN Refer to Section **2.2.3A**.

SYNC OUT Refer to Section **2.2.3A**.

2.2.4 POWER REQUIREMENTS **+24V @ 30 mA** **+12V @ 125 mA**
-24V @ 0 mA **-12V @ 0 mA**

2.2.5 OTHER INFORMATION

WEIGHT (SHIPPING) 3.2 lbs. (1.45 kg)
(NET) 2.1 lbs. (.95 kg)

DIWENSIONS Standard single width **NIM** module (1.35 x 8.714 in.) per **TID 20893** (Rev.).

INSTRUCTION MANUAL One-provided with each instrument.

3.0 INSTALLATION

3.1 TC 178

3.1.1 POWER

The TC 178 Quad **Preamplifier** is not self-powered and must be connected via the power cable to a main amplifier with provisions for providing preamplifier power or a separate preamplifier **power** supply. Refer to the CAUTION at the **beginning** of the manual before connecting the TC 178 power cable to TENNELEC main amplifiers other than the TC 240 series.

3.1.2 CONNECTIONS

3.1.2A DETECTOR CONNECTION

To preserve the low-noise characteristics of the system, the capacitance to ground at the input of the preamplifier should be **kept** to a minimum. If a cable between detector and preamplifier must be used, it should be as short as **possible and** it must be shielded with one end of the shield connected to the detector housing and the other to the preamplifier housing. Double shielded cable (RG71/U) is preferable to single shielded cable (RG62/U).

To avoid microphonics, it is desirable that the geometrical relationship between detector and preamplifier be kept rigidly fixed.

The length of cable connecting the preamplifier and the detector should be kept to a minimum for reasons of stability in addition to noise **considerations**. The cable connecting the **preamplifier** and detector introduces a phase shift into the preamplifier **feedback** loop which adversely affects stability. A maximum length of cable cannot be assigned to the large number of **detector** and cable **combinations**, but a typical maximum length for the TC 178 is 10 ft.

The noise performance of the preamplifier can be estimated from the sum of connecting cable and detector capacitance. The noise as a function of this input capacitance is shown in Figure 0 2.1 and 2.2 for peaking times of 4 **usec** and 1 **usec**.

3.1.2B BIAS SUPPLY CONNECTION

In the TC 178, the bias connection is routed through the preamplifier case.

If a battery pack is used for bias, **no** special precautions need be **taken in** cable routing to avoid noise pickup provided the case of the battery pack is connected to **the preamplifier**.

If a power line operated supply is used, or if a battery pack is used which is grounded to the main amplifier frame, then it **is** desirable to take the following precautions to avoid ground-loop pickup.

- a. Locate the power supply physically close to the main amplifier.
- b. Ground the supply to the **main** amplifier with large-gauge wire or shield braid at least **1/4" (6mm)** wide.
- c. Cut the high voltage cable to approximately the same length as the preamplifier signal cable and twist or taps the two together.
- d. Never plug the main amplifier and HV supply into different wall outlets. If necessary, use a local **distribution** box for all components of a **spectrometer** system to avoid **making** the building part of a ground loop.

3.1.2C ENERGY OUTPUT CONNECTION

The energy output (**E OUT**) **is intended** to drive a 50 ohm line (**RG 58/U**) **which** may be connected directly to the input of the main amplifier. A 50 ohm termination **is** not required **as the preamplifier** is stable unterminated. The preamplifier will drive any length of cable; **however**, for long cable lengths, cable losses must be considered..

To minimize ground-loop noise pickup, the following pattern of connections **is recommended**.

- a. Place **the test** pulse generator as close to the main **ampifier** as possible.

- b. Cut the signal and test pulse cables to approximately the same length as the $\pm 12V$, $\pm 24V$ wires. Twist or tape all cables into one bundle.

The purpose of instruction.(a) is to prevent power supply noise spikes (which nearly always exist between widely spaced ground points **in a NIM** bin) from appearing in series with the signal ground returns. The purpose of instruction (b) is to avoid local radio station pickup which frequently occurs because of the loop-antenna effect in a network of spread-out cables.

The TC 178 **E OUT** signal is terminated at the sending end by resistor **R105** (CHANNEL A) R129 (CHANNEL B), R131 (CHANNEL **C**) and R133 (CHANNEL **D**) as shown on the TC 178 schematic. In the TC 178, these resistors are nominally 49.9 ohms each to match 50 ohm cables. Other cable impedances may be used by changing to resistors of appropriate values.

3.1.2D TEST IN CONNECTION

The TC 178 TEST IN connector (one per channel) may be used in conjunction with a TC **178P Quad** Pulser to operate the TC 178 internal pulser. The TC 178P controls the **rate and** amplitude, the polarity is controlled by the TC 178 PULSER POLARITY control (one per channel). The TC 178 TEST IN **connector** CANNOT be used in conjunction with a standard pulse generator. The signal to the TEST IN connector must provide a trigger pulse (TTL level) and a dc level to provide the pulse amplitude information.

If the TC 178 is used without a TC **178P**, a test pulse can be injected at the TC 178 INPUT by using an external test capacitor and terminating resistor as shown in Figure 3.1. The use of an external test capacitor will give a more accurate representation of the TC 178 transient response than using the TC 178 internal pulser.




Figure 3.1 External Test Input Terminator

3.1.3 GENERAL PRECAUTIONS

3.1.3A FET PROTECTION

In the TC 178, a FET connected as a diode and a resistor in series with the input are used to protect the input FET against accidental short circuits. The protection network **degrades** the preamplifier noise slightly (refer **to Section 5.1** and **Figures 5.1 through 5.3** for **more detail**). The **major noise source** of the protection network (especially at short peaking times, and detector capacitance of **100pF** or greater) is **resistor R4** (47 ohm). This resistor is easily removed or replaced with a lower value to retain some protection with improved preamplifier noise performance.

IMPORTANT: See Section 7.1 for instructions on how to properly open or remove the TC 178 case.

3.1.3B APPLYING BIAS VOLTAGE

In the following **statements, it is** essential that the user recognize the **distinction between** rapid voltage changes at the B.V. IN connector and the signal INPUT connector.

The TENNELEC TC 178 preamplifier can safely withstand the application of detector bias voltage in **±500V** steps spaced **10 sec** apart, with or without protection for the FET, if **the** voltage is applied through the H.V. IN connector.

Without protection, short circuits at the signal INPUT terminal, may cause FET damage if the bias exceeds **50V**. **Connecting a preamplifier** to a detector with bias **voltage** applied, either through the preamplifier or directly to the detector, is nearly equivalent to a short circuit.

With the protection **network** installed, an occasional short circuit at the signal INPUT terminal will not **cause** FET damage if the bias is 500V or less.

With or without protection, the preamplifier may be disconnected **from a** charged detector, but not reconnected except within the **limits** stated above for short-circuiting the **signal** INPUT terminal. When a preamplifier is disconnected from a charged detector,.

the signal; **INPUT** terminal should not be short-circuited except **within** the limits **given** above. Additionally, **the detector bias supply should not be disconnected without first** reducing the voltage to zero in 500V steps spaced 10 **sec** apart and then waiting for an additional minute **to allow** the preamplifier filter **network** to discharge. **The reason** for this last precaution is that **without** a **return** path through the **power supply, neither the filter capacitors nor the input coupling capacitor will have a discharge path; it may take an hour for discharge to occur through leakage resistance alone.**

The **user is** reminded at this point that because of their vulnerability to accidental damage, the **FETs** are not covered by warranty.

3.1.3C DETECTOR BIAS NETWORK

The TC 178 detector **bias** network **consists** of two 1 **megohm** resistors in the filter network and one 1 **megohm** resistor for, **the detector load resistor.**

If the detector load **resistor** and bias resistors supplied as **standard** are **inappropriate**, the user can install **different** values. The noise performance of the **preamplifiers** will **not** improve significantly if the detector load resistor is **increased** unless the **charge-loop feedback resistor is increased.** For **modifications** of this nature, see Section 7.3.

3.1.3D DETECTOR BIAS LEVEL

If a system containing a **room-temperature surface-barrier** detector is assembled and turned on in the absence of bias voltage, the **electronic** noise will be high. As **the bias** voltage is **increased**, the noise (observed as "grass" on an oscilloscope connected to the output of the main amplifier) should drop sharply. As the voltage is increased further, the noise should continue to drop up to the point where rated detector voltage is reached, then it should increase again. The appearance of this **noise** at this higher than rated-bias level will **not be clean "grass"** as observed earlier, but will appear as a **series of discontinuities** on the baseline. This later appearance is characteristic of avalanche **breakdown** in the **detector.** The correct operating voltage for the **detector** is about 10% **below** this avalanche level. If the load resistor is too high,

for the leakage current of the detector, the noise level, will not drop as the bias voltage is increased until it exceeds the IR drop in the load resistor. In extreme cases, this will not occur until the bias voltage as indicated on the power supply exceeds the detector manufacturer's specified maximum. **NOTE:** To compensate for the voltage drop in the detector bias network, the total dc resistance is 3 megohm. This resistance times the detector bias current will give the voltage drop across the detector bias network. The voltage applied to the detector is the power supply voltage minus the voltage drop across the detector bias network.

The detector voltage is given by

$$V_{DET} = V_{BIAS} - (I_{LEAKAGE} \times R_{LOAD})$$

where V_{BIAS} is the output voltage to the high voltage power supply, $I_{LEAKAGE}$ is the detector bias current and R_{LOAD} is the total resistance in the preamplifier detector bias network (3 megohm for the TC 178).

CAUTION: The user is encouraged to discuss the detector break-down characteristics with the manufacturer of the detector. TENNELEC cannot assume responsibility for damage to the detector caused by improper load-resistor selection or by improper application of detector bias voltage; With some detectors, permanent damage will result from over-voltage.

With time and radiation damage, the onset of avalanche noise (also known as "flicker" noise), may drop to a level below the ratings of the detector, requiring a reduction of operating voltage for acceptable energy resolution (or background count rate).

3.2 TC 178P

3.2.1 POWER

The TC 178P QUAD PULSER requires a standard NIM bin and power supply, such as the TB3/TC 911, for operation. The bin provides mechanical mounting and power supply distribution. Always turn OFF the bin power supply when inserting or removing any modules.

The TC 178P is designed so that it is not possible to overload the power supply, even with a full compliment of modules in the bin. **Since** this may not be true when the bin contains modules other, **than those** of TENNELEC design, the **power** supply voltages should be checked after all modules **have** been inserted. The **TENNELEC** Bin and Power Supply provides power supply test points on the bin control panel for monitoring the dc voltage levels.

The TC 178P PULSE AMPLITUDE accuracy does not depend upon the **+24V** power supply level or level changes except to a very 'minor degree (**<.1 mV** per Volt variations of the **+24V** level).

The accuracy of the TC **178P** RATE control is dependent on the **+24V** power supply level. The TC 178 RATE varies approximately directly with changes of the **+24V** level. For optimum RATE accuracy, the **+24V** power supply level should be **+23.9 to +24.1V**. The **TB3/TC** 911 Nim Bin and Power Supply is **furnished** with trimming controls for precisely setting the voltages., On TENNELEC bin supplies, all bias voltages are of instrument-grade quality, highly regulated and very stable.

3.2.2 CONNECTIONS

3.2.2A PULSER OUT

Each channel of the TC 178P has a PULSER OUT connector on the rear panel. The PULSER OUT signal 'is intended to be connected to the TC 178 **TEST** IN connector via coax cable such as **RG58/U**. The PULSER OUT signal consists of a TTL signal superimposed on a dc level. The PULSER OUT signal **is** separated into its components in the TC 178.

DO **NOT TERMINATE** THE **PULSER OUT** SIGNAL. Coaxial cables other than **RG58/U** may be used; however, the cable characteristic impedance should be 50 ohms as the TC 178P PULSER OUT **signal** is series terminated in **50** ohms.

Do not connect the TC 178P PULSER OUT to the TC 178 Preamplifier **using a** test capacitor. The TC 178P is not designed to be used in this manner.

3.2.2B SYNC OUT

The TC 178P has a SYNC OUT connector on the front panel and the rear panel. The outputs are identical and driven from the same point; The output signal (see Section 2, Figure 2.5 for a typical waveform) is a 60 nsec **wide** positive pulse. The output is series terminated and the output impedance is 50 ohms. The SYNC OUT signal can be used to trigger oscilloscopes or gate instruments to allow observation of the test signal from the TC 178. When using the SYNC OUT signal as a trigger or gate signal, the signal should be connected using 50 ohm characteristic impedance coaxial cable such as **RG58/U** or **RG174/U** and be terminated into 50 ohms.

3.2.2C EXT IN

The TC 178P has an EXT IN connector on the front panel and the rear panel. The inputs are identical and drive the **same** point. The **input** signal requirements are **≥ +2V, ≥100** nsec wide. A typical EXT IN signal would be a TTL or NIM logic pulse, approximately **+5V** and 100 nsec to 500 nsec wide. The EXT IN signal can be used to trigger the TC **178P** from a random source or to synchronize several TC **178Ps** by use of an external generator or synch'ronize several TC **178Ps** to a master TC 178P.

4.0 OPERATING PROCEDURES

4.1 FIRST-TIME OPERATION

Users will find **it helpful** to familiarize themselves with the TC 178 QUAD' PREAMPLIFIER and TC **178P** QUAD PULSER by conducting a few simple tests.

4.1.1 EQUIPMENT REQUIRED

1. NIM Bin and Power **Supply** (**TENNELEC TB3/TC 911** or equivalent).
2. Oscilloscope (**TEKTRONIX 465** or equivalent).
3. Shielded 50 ohm cables with **BNC** connectors;
4. **BNC TEES** and 50 ohm terminators.
5. Vacuum chamber (**TBNNELEC TV-1** or equivalent).

6. Charged Particle detector (TENNELEC Surface Barrier Detector, or equivalent, size not critical).
7. Alpha **source** such as **Am²⁴¹** or **Po²¹⁰**.
8. High Voltage Bias Supply (TENNELEC TC 908 or equivalent).
9. Shaping Amplifier (TENNELEC TC 243 or equivalent).
10. Multichannel Analyzer (NUCLEUS Spectrum **88** or equivalent).

4.1.2 INITIAL TEST SETUP

Set the TC 178 controls as follows:

RANGE	0.1 GeV
RISETIME	FULLY COUNTERCLOCKWISE (use insulated screwdriver)
DC OFFSET	NO ADJUSTMENT REQUIRED
PULSER POLARITY	NEG

Set the TC 178P controls as follows:

RATE	1.0k
MULT	10.0
PULSE AMPLITUDE (CHANNEL A - D)	10.0
ON/OFF (CHANNEL A - D)	ON

Set the **OSCILLOSCOPE** controls as follows:

VERTICAL MODE	CH1
CH1 VERT SENS	0.5V/DIV (dc-coupled)
HORIZ SWEEP	50 nsec/DIV
TRIGGER	AC/NORM/EXT/+

The TC 243 is used **as a source of preamplifier** power in this test **and the control settings are** not important.

Connect the TC 178 **power cable** to the TC 243 preamplifier **power connector**.

Install the TC 178P and TC 243 in the NIM Bin and Power supply and apply power.

Connect **the** equipment as shown in Figure 4.1.




Figure ~4.1 Initial Test Setup

Observe a waveform on the OSCILLOSCOPE similar to Figure 4.2. Vary the **RISETIME** control (use an insulated blade screwdriver) for best transient response. The amplitude of the waveform should be **+1.9** to **+2.1V**. If the **amplitude** is outside these limits, use the 0.1 **GeV** PULSER CALIBRATION control to obtain a **+2.0V** amplitude (use an **insulated** blade screwdriver).




Figure 4.2 Typical TC **178** Output, 0.1 Gev Range

Set the TC 178 **RANGE** switch to 1.0 GeV and adjust the **RISETIME** control (use an **insulated** blade screwdriver) for best transient response. The waveform should be similar to **Figure 4.3** and the amplitude should be **+1.9 to +2.1V**. If the amplitude is outside these limits, use the 1.0 GeV **PULSER CALIBRATION** control to obtain a **+2.0V** amplitude (use an **insulated** blade screwdriver).




Figure 4.3 Typical TC 178 Output, 1.0 GeV Range

Set the TC 178 **RANGE** switch to 10 GeV and adjust the **RISETIME** control (use an **insulated** blade screwdriver) for best transient response. The waveform should be similar to **Figure 4.4** and the amplitude should be **+1.9 to +2.1V**. If the amplitude is **outside** these limits, use the 10 GeV **PULSER CALIBRATION** control to obtain a **+2.0V** amplitude (use an **insulated** blade screwdriver).



Figure 4.4 Typical TC 178 Output, 10.0 GeV Range

The RISETIME, control has greater range, and is more sensitive for the 1.0 GeV and 10 GeV energy ranges. The risetime range can be verified for a detector and cable combination by installing a coaxial cable of the same length used to connect the detector and preamplifier and a capacitor that represents the equivalent detector capacitance on the TC 178 INPUT. This will allow easy evaluation of any problems associated with cable length between the detector and preamplifier.

By decreasing the TC 178P PULSE AMPLITUDE control, the trigger point of the TC 178 E OUT LEDs can be verified. This level should occur at approximately 2% of full scale. Full scale is 0.1, GeV, 1.0 GeV and 10 GeV and the energy of the test pulse is indicated by the PULSE AMPLITUDE dial.

The response of the TC 178 is optimized for output signals in the positive direction. However, operation in the negative direction is easily verified by changing the TC 178 PULSER POLARITY to positive (PGS). The internal pulser is optimized for use in the NEG position (positive preamplifier output). The negative transient response of the TC 178, using the internal pulser in the POS position, will not be as accurate as the TC 178 positive transient response. To properly evaluate the transient response in either direction, an external pulser and test capacitor should be connected to the TC 178 INPUT.

NOTE: The internal pulser of the TC 178 provides a rectangular pulse. The initial edge is used as the input pulse. The trailing edge of the rectangular pulse will produce an output pulse of the opposite polarity of the initial output pulse. The trailing edge pulse will occur approximately 40 usec after the leading edge.

4.1.3 PULSBR CALIBRATION

To properly calibrate the TC 178 pulser requires a detector and radioactive with either a primary peak such as Am^{241} or a monoenergetic source such as Po^{210} . The use of natural alpha sources will only allow calibration on the lower 10% of the 0.1 GeV range.

When the TC **178** and TC 178P are used at higher energy8 a8 with an accelerator, the ability **to calibrate** to a beam energy will allow greater flexibility in verifying overall system calibration.

4.2 **ADJUSTMENTS**

4.2.1 **TC 178**

4.2.1A **DC OFFSET VOLTAGE**

Each **channel of** the TC 17.8 has an adjustment for the dc **offset** voltage as **measured** at the TC 178 E **OUT**. The DC OFFSET adjustment of the TC 178 is accessible through a **hole (one per channel)** in the preamplifier case. This hole is normally covered with a press-fit **plug** to reduce electrical pickup. The DC OFFSET potentiometer is R23 a8 shown **on the** TC 178 schematic included in **this manual**, **This adjustment** affect8 the drain current of the input FET (**Q1**) and therefore the noise performance.

To ensure the optimum **noise performance**, the DC offset should be adjusted for **-100 mV** as measured at the energy output connector. This measurement is made with the energy output connector unterminated. The DC offset voltage of the TC 178 is adjusted **to -100 mV (±50 mV)** before leaving the factory and should not require further adjustment unless the input **FET** is replaced.

4.2.1B **RISETIME**

Each channel of the **TC 178** has a **risetime** adjustment. ALWAYS USE AW **INSULATED** SCREWDRIVER (either plastic shaft or covered metal **shaft**). The **risetime** adjustment correspond8 to **R21 on** the TC 178 schematic. To obtain optimum transient response and timing performance, the **risetime** adjustment must be set with the detector (or an equivalent detector capacitance) and any connecting cable connected **to the** preamplifier input and bias applied (if an equivalent detector capacitance is used in place of the **detector**, no high voltage bias is **necessary**). When **minimum risetime** is obtained, a short period overshoot of 10% **to** 20% will be observed. This overshoot ha8 **no effect** on noise performance or linearity within the normal dynamic range.

Additionally, the overshoot will have no adverse effect on timing measurements. **Further** improvement in timing performance can be obtained by adjusting the **risetime** control for a shorter **risetime** and more overshoot. The increased overshoot will be accompanied, by some ringing which limits the minimum risetime. Some **ringing** can be tolerated if the discriminator has an adjustable **deadtime** control.

NOTE: Detector inductance has a very important bearing on preamplifier signal **risetime** and **circuit** stability. It is essential when operating the preamplifier for the first time to match it to the specific detector characteristics. To do this, connect a fast oscilloscope to the output of the preamplifier and excite the detector with a weak **radiation source** that generates only one spectral line. **Po²¹⁰** is a good source for this purpose with room temperature silicon detectors. The oscilloscope should have a bandwidth of 250 **MHz** or more, and should be triggered internally. A **continuously** changing **ensemble** of randomly occurring pulses will be seen, making observation of the pulse difficult. Unfortunately, not **much** can be done about this.

The adjustment is made by turning the **risetime** control to obtain the optimum balance between ringing and risetime.

Depending on the equivalent circuit of the detector, some settings may cause oscillation. This is normal. A stubborn case of oscillation indicates excessive detector inductance or connecting cable length. If the oscillation is difficult to eliminate, or elimination results in slower than expected risetime, it is recommended that the connecting **cable** be reduced to as short as practical. The **TC 178** normally will be stable with up to 10 feet of cable between detector and preamplifier.

Typical **risetime** versus detector capacitance is shown in Section 2.1. Due to the fast risetimes of the **TC 178**, it is **recommended** that a **pulsar** and oscilloscope with 1 nsec or less **risetime** be used to check and adjust the **transient response** when not using a detector and source. The actual preamplifier **risetime** (t_{rpa}) is

$$t_{rpa} = (t_{ra}^2 - t_{rb}^2)^{1/2}$$

where t_{ra} is the **risetime** of the **preamplifier, pulser** and **oscilloscope**, and t_{rb} is the **risetime** of the pulser and **oscilloscope** only., By using the equation given above, the **risetime** contributions of **the pulser** and **oscilloscope** can be removed. All **risetime** specifications apply when the output is terminated in its characteristic impedance (50 ohms).

4.2.1C PULSER CALIBRATION

The TC 178 QUAD PREAMPLIFIER is shipped with the PULSER CALIBRATION controls set to deliver nominally **+2.0V** E OUT pulses for a **10.0V** (Full Scale) setting on the TC 178P PULSE AMPLITUDE control (TC 178 PULSE POLARITY set to NEG.). This setting (in terms of energy) should be within **±5%** for the selected energy range. The largest error will be for the 0.1 **GeV** range as stray feedback capacitance and feedback capacitor tolerance represent a larger percentage error for this range. The TC 178 **PULSER CALIBRATION** can be accomplished by either of two methods. One method requires a source of known energy and a detector., The energy should be greater than 10% of the full scale range to allow for greater accuracy in calibration (minimum energy of 10 **MeV**, 100 **MeV** and 1 **GeV**). The second, and more practical, **requires** injecting a known amount of charge into the input of the TC 178 and **adjusting the TC 178 PULSER CALIBRATION** for a corresponding output amplitude. The actual preamplifier output voltage or a shaped signal can be used for the **adjustment**. For instance, given a **1pF** test capacitor, to inject a charge equivalent to 10 **MeV (Si)** requires a **0.439V** signal, 100 **MeV (Si)** requires a **4.39V** signal. To **generate** a 1.0 **GeV** signal, a test capacitor of **10pF** requires a **4.39V** signal. **The 10 GeV** signal can be simulated by using a **100pF** test capacitor and a signal voltage' of **4.39V**. Use of test capacitors **over 100pF** are not recommended. The smaller the test capacitor, the less are the effects on transient response and noise. The use of a nonmetallic screwdriver is recommended **for adjusting the PULSER CALIBRATION** control. The controls are accessible from the bottom of the TC 178 and the proper sequence of adjustment is 0.1 **GeV**, 1.0 **GeV**, and 10 **GeV**. There is **some** interaction and the adjustment sequence should be repeated until sufficient accuracy is obtained.

4.2.2 Tc **178P**

No user adjustments, other than front panel controls, are required.

5.0 NOISE

NOISE PERFORMANCE

To convert from full-width-at-half-maximum (**FWHM**) Si to a different reference, use Table 1.

TABLE 1

CONVERSION OF FWHM Si TO OTHER REFERENCE VALUES

Reference	Multiply eV FWHM Si by:
FWHM Si	1.00 (W=3.6 eV/electron-hole pair)
FWHM Ge	0.819 (W=2.95 eV/electron-hole pair)
FWHM P10	6.94 (W=25 eV/electron-hole pair)
Ion pairs rms	0.144
Coulombs rms	2.3×10^{-20}

The noise performance of the TC 178 for a peaking time of 1.0 **usec** and 4.0 **usec** is given in **Section 2.1**. The noise performance of the TC 178 at various shaping times and detector load resistance, with and without **the** protection network, is shown in **Figures 5.1** through **5.3**. With the aid of this data, the noise performance of the preamplifier can **be predicted** for almost any combination of detector capacitance and shaping time. It is stressed that the foregoing figures are noise levels and not spectral resolution. The final spectral resolution depends not only on the preamplifier noise but also on the type of detector used, the count rate, and other factors. An additional consideration in evaluating preamplifier noise limitations is the detector leakage current. All the previous noise data is **representative** of the preamplifier detector combination with zero detector leakage current.




Figure 5.1 Electronic Noise vs Equivalent Detector Capacitance




Figure 5.2 Electronic Noise vs Equivalent Detector Capacitance




Figure 5.3 Electronic Noise vs Equivalent Detector Capacitance

5.2 NOISE MEASUREMENTS

To verify the proper operation of the preamplifier, *noise* measurements can be made by either of two methods. These **measurements** can be compared with the values given in Section 2.1 and Figure 5.1 through 5.3.

One method requires a calibrated step-generator (e.g. TENNELECTC812), a shaping main-amplifier (TC 205A, TC 240, TC 241, TC 242, TC 243, or TC 244), and a multichannel pulse height analyzer. After the shaping time constant has been chosen and the analyzer has been calibrated in terms of energy per channel, pulses are fed through the test capacitor and the line width recorded by the analyzer is measured. For this test, as for any measurement of absolute noise of the preamplifier, the detector **should be replaced** by a **dummy** capacitor of the **same** capacitance. The FWHM of the line should be close to the values **given** in Section 2.1 for typical performance and the typical data given in Figures 5.1 through 5.3. If the noise at 4.0 **usec** peaking time exceeds the guaranteed values given in Section 2.1, verify that the correct detector load resistor is installed.

The second method requires the use of a calibrated pulse generator, a shaping main-amplifier, an **average-type** ac voltmeter (such as a Hewlett-Packard 400D, 400H, or 400L) or a true rms voltmeter (such as a Hewlett-Packard 3400A), and a calibrated oscilloscope. A step of known **amplitude** V_i is applied to the input through the test **capacitor** C_T , resulting in a charge transfer to the input of $V_i \times C_T$ coulombs. The resulting main-amplifier **pulse height** V_o is recorded with the oscilloscope. The **pulse generator** is then turned off, **and** the true rms noise level V_n is **measured** at the output of the main amplifier; **if** a true rms voltmeter is used, the reading is directly V_n . If an average-type voltmeter is used, the reading V_r should be multiplied by the **factor** 1.128 to obtain V_n . The level in **keV** FWHM referred to Si detectors **is given** by

$$\text{Noise (FWHM)} = \frac{V_n \times V_i \times C_T}{V_o} \times 5.298 \times 10^{16}$$

where 5.298×10^{16} is a factor that contains the charge of an electron in coulombs, the **energy** necessary to produce one electron-hole pair in **silicon**, and the conversion constant between rme and **FWHM**.

For detectors, other than silicon, choose the appropriate multiplier, from Table 1.

If a problem with excessive noise should occur, either of the two procedures described should be used to evaluate the noise performance. The preamplifier noise performance can be verified by replacing the detector with a suitable capacitor having the same capacitance value as the detector. If this noise is within specifications at 4.0 usec peaking time, or similar to data given in section 5.1, the problem is associated with the detector. The total noise of the system is given by

$$N_{\text{total}} = \{(N_{\text{preamplifier}})^2 + (N_{\text{detector}})^2\}^{1/2}.$$

Using the above equation and the noise of the preamplifier (as previously determined) the noise of the detector can be calculated and compared with the manufacturer's data.

6.0 COUNT RATE EFFECTS

6.1 RESOLUTION

The shape of a typical spectral line is Gaussian as shown in Figure 6.1.

The resolution, or ability of a nuclear spectrometer to separate different radiation energies, is usually expressed in terms of the full width of the spectral lines measured at half their maximum height. This quantity is denoted by the letters FWHM and is given in units of energy. The FWHM is 2.35σ where σ is the standard deviation. If noise alone controls the resolution, then σ and the rms noise level are synonymous. We have been using FWHM to characterize preamplifier noise levels in the previous section.




Figure 6.1. The shape of a typical spectral line (Gaussian).

The resolution obtained in any particular spectral measurement is the result of several factors: preamplifier noise, detector characteristics, count rate, radiation energy, overall system stability, proper interfacing between instruments within the system, etc. In an experimental situation in which the count rate is low enough so that the pulse shape can be adjusted for the best signal-to-noise ratio without being affected by pile-up or -baseline shift but high enough so that effects due to long term drifts can be neglected, the resolution will be determined by three factors: (a) detector resolution for the particular radiation energy being observed, (b) electronic noise, and (c) interfacing. Furthermore, if it is assumed that the different components of the system are properly matched, the line-width is a function of only the detector resolution and the electronic noise. The two noise components are related in the following way:

$$R^2 = (\text{Total Resolution})^2 = (\text{Detector Resolution})^2 + (\text{Electronic Noise})^2$$

We shall call R the intrinsic resolution of the system. In a counting situation in which the conditions are not ideal, the measured resolution will be worse than the intrinsic resolution. Usually, the main factor in line-width broadening is count rate. Count rate can have a deleterious effect in spectral resolution through several mechanisms. The three most commonly found are pile-up of pulses, baseline shifts, and thermal effects in components. The last two can usually be neglected in properly designed systems; the first one is more difficult to contend with.

Usually, pile-up of shaped pulses in the main amplifier will set the practical upper count rate limit. However, at very high energy (lowest gain settings of the main amplifier), the limitation may occur in the preamplifier. A discussion of preamplifier pile-up follows, plus the technique of computing the upper count rate limit. The pulse obtained at the output of the preamplifier appears as shown in Figure 6.2.




Figure 6.2 Pulse Shape at the Output of a Preamplifier

When pulses come in rapid succession, the wave form at the output of the preamplifier appears as shown in Figure 6.3. The dotted line at the top of Figure 6.3 indicates the limit of the linear range of the preamplifier.

If the count rate is high enough, some of the pulses will rise beyond the linear range and therefore, their amplitudes will be distorted. (The meaning of "linearity" is explained in Section 6.2.) If we assume 10V to be the limit of the linear range, the preamplifier sensitivity to be 20 mV/MeV, the average radiation energy to be 100 MeV, and the decay time constant to be 50 usec, we can compute the count rate that will be necessary to make 5% of the pulses fall beyond the linear range from the formula

$$\text{where } n = \frac{2}{T_d} \left[\frac{V_m - E G_c}{2.5 E G_c} \right]^2$$

n = count rate in cps .

T_d = decay time constant in sec.

V_m = linear range in volts

E = radiation energy in MeV

G_c = preamplifier sensitivity in V/MeV.

Replacing symbols by actual numbers,

$$n = 102 \times 10^3 \text{ cps @ } 100 \text{ MeV}$$

Since the TC 178 charge loop is ac-coupled, a count-rate product cannot be assigned with any useful units as the number would apply only for one specific energy. The minimum count rate at full scale for each range of the TC 178 is given in Section 2.1.




Figure 6.3 The Waveform at the Output of a Preamplifier with Pulser Applied to the Input in Rapid Succession

6.2 NONLINEARITY

If a graph of output, signal_level_V , vs. input pulse height: V_i (dynamic **characteristic**) is drawn, a perfectly straight line passing through the origin should result. In practice, the dynamic characteristic could have a slight curvature up to a certain signal level, **beyond** which the curvature increases **drastically**. The onset of this drastic change is usually considered to be the upper limit of the normal dynamic range (rated. output).

Integral nonlinearity is defined **as the** maximum **deviation** of the **measured preamplifier** response from the ideal response, expressed as a **percentage** of the rated output (as described in **the** preceding. paragraph). This definition is useful only for isolated preamplifier pulses as shown in Figure 6.2. When pileup occurs as the result of an ensemble of closely spaced small pulses (Figure 6.3), we are interested not only in the integral nonlinearity **but also** in the deviation of height of **individual** steps (within the 'linear' range **of the preamplifier**) from the expected height. This incremental deviation in $\Delta V_o / \Delta V_i$ from the value at zero volts on **the dynamic** characteristic is described as the **differential nonlinearity**. It is this definition which is **used in** the table of specifications.

7 . 0 PREAMPLIFIER **MODIFICATIONS**

REMOVING THE CASE

Remove the 24 press fit plugs, **the** four mounting screws on the bottom of **the case**, and carefully **remove** the preamp from the case.

7.1 CHANGING PREAMPLIFIER SENSITIVITY

Reducing the preamplifier sensitivity will almost certainly cause it to oscillate unless the stabilization networks are changed as well. For this change the user is requested to return the instrument, to TENNELEC for modification. **In**creasing the sensitivity will not cause oscillation, but may degrade the **risetime** and pulse shape,

7.2 FRONT-END MODIFICATIONS

The TC 178 is supplied with one 1 megohm resistor for the **detector** load resistor and two 1 megohm resistors in the detector **bias** voltage filter **network**. If these values **are** not appropriate (too low), the 1 megohm detector. **load** resistor should be increased. Do not increase the detector bias voltage filter network resistors.

If any resistors in the detector bias network are changed, it is very important that all solder joints in the high voltage chain be smooth-surfaced and have no protruding sharp points. Furthermore, it is important that all capacitors and high megohm resistors used in this part of the **circuit** be free of surface **contamination**. Components **that** are contaminated can cause increases in preamplifier noise, leakage current, noise spikes from arcing, etc.

7.3 MODIFICATION OF FET DRAIN CURRENT

The FET drain current of the TC 174 can be changed by adjusting **R8**. **However**, this will affect the dc offset, which should never be set to less than -50 mV. The normal dc offset voltage is -100 mV, which, is close to the value resulting from operating the input **FET** near **I_{dss}**. Best noise **performance** is obtained at or near the **I_{dss}** level.

7.4 ELECTRONIC NOISE REDUCTION

The noise performance of the TC 178 can be changed significantly by eliminating the series **protection** resistor, increasing the **detector** load to 10 megohm, and changing the feedback resistor to a higher value such as 100 megohm. Changes of this nature are best performed at TENNELEC.

8.0 SHIPPING DAMAGE

Upon receipt of the instrument, examine it for shipping damage. **Damage claims should** be filed with the carrier. The claims **agent** should **receive a full report**: a copy of that report should be sent to TENNELEC, Inc., P.O. Box **2560, Oak Ridge, Tennessee** 37830-2560. The model number and serial number of the instrument must be included in the report. **Any** remedial action taken by TENNELEC, Inc., will be based on the information contained in this report.

9.0 SERVICING

In the event of a component failure, replacement may be done in the field or the instrument may be returned to our plant **for repair**. There will be no charge for repairs that fall within **the warranty**.

10.0 WARRANTY

In connection with **TENNELEC's** warranty (inside front cover), **TENNELEC** suggests that **if a** fault develops, the customer should immediately notify the **TENNELEC Customer Service Manager**. He may be **able** to prescribe repairs **and send replacement parts** which will enable you to get the instrument operating sooner and at less expense than if you returned it.

Should return prove necessary, the **TENNELEC** Customer Service Manager must be informed **in WRITING, BY CABLE** or TWX of the nature of the fault and the model number and serial number of **the** Instrument. Pack the instrument, well and ship **PREPAID** and **'INSURED** to **TENNELEC, Inc., 601 Oak Ridge Turnpike, Oak Ridge, Tennessee 37830-2560**. As stated in the warranty **DAMAGE IN TRANSIT WILL BE REPAIRED AT THE SENDER'S EXPENSE** as will damage **that obviously** resulted from abuse or misuse of the instrument.


Quotations **for** repair of such damage will be sent for your approval before repair is undertaken.

* * * * *
*
* **TENNELEC'S** Quality Assurance Program **re-** *
* **quires** that each and every instrument be *
* **fully** aged, vibrated, and electronically *
* **checked.** *
*
* Should the user require a copy of the *
* **Quality** control Procedure and Test Record, *
* **please** call the Customer Service **Depart-** *
* **ment** of TENNELEC. Both model number and *
* **Serial** number are required. *
*
* * * * *

NANUAL REF.: 0

3/85 - Engineering and component improvements may be made after date of printing.

REV	ECN	REVISIONS	DATE
0			8-27-59




TÉNNÉLEC
 A DIVISION OF
 GENERAL ELECTRIC

TC 178P
 CIRCUIT DIAGRAM
 PRELIMINARY

DATE: 8-27-59
 BY: J. S. W.
 CHECKED: J. S. W.
 APPROVED: J. S. W.
 TC 178P

CHANNEL A



PRELIMINARY

TÉNNELÉC
DATA MODEL: TTNM178A

TC 178
CIRCUIT DIAGRAM

REV	DESCRIPTION	REV	DATE
1	5.3.84/8.20/11.1	1	11.1.84
2	8.11.84/8.20/11.1	2	11.1.84

