

INSTRUCTION MANUAL

TC 248 FAST/SLOW AMPLIFIER

WARRANTY

TENNELEC, INC. warrants that the products or components manufactured by it shall be free **from** defects in material or workmanship **for** a period of two years from the date of delivery to purchaser. If such product or component is determined to be defective by TENNELEC, its sole warranty obligation shall be limited to either replacing or repairing such defective product or component or allowing credit therefor, at TENNELEC's option. Such warranty is further conditioned upon the purchaser's giving prompt notice of any such defect and satisfactory **proof** thereof to TENNELEC's customer service manager, thereafter upon TENNELEC's approval, the purchaser shall return such defective product or component to TENNELEC's factory at Oak Ridge, Tennessee, all transportation charges prepaid. TENNELEC shall be responsible only for transportation charges incurred in returning **such** product **or** component to purchaser. All customs, brokerage and duty charges shall be at the expense of the purchaser. Damage in transit due to inadequate packaging will be repaired at purchaser's expense. Any repairs or replacements by the purchaser without TENNELEC's approval, any willful abuse or any evidence that the product or component was not properly used and maintained, would automatically void this warranty.

TENNELEC makes no warranty whatsoever in respect to products or components not manufactured by it but instead the applicable warranties, if any, of the respective manufacturers thereof shall **apply**. Likewise fuses, batteries and input transistors in ultra low-noise amplifiers are specifically excluded from this warranty.

THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, STATUTORY OR OTHERWISE, INCLUDING WARRANTY OF MERCHANTABILITY AND FITNESS.

*** * * * * CAUTION * * * * ***

The **PREAMP** connector on this instrument is directly compatible only with TENNELEC preamplifiers with serial numbers greater than 2000. It is also directly compatible with standard Aptic, Canberra, EG&G Ortec and PGT preamplifiers.

If a TENNELEC preamplifier with serial number less than 2000 is used, then a model ADT 1 PREAWP POWER ADAPTER must be used.

If there are any questions regarding the compatibility of the **PREAMP** connector of this instrument, please contact the TENNELEC Marketing Department for assistance.

*** * * * * WARNING * * * * ***

*** Improper connection to the **PREAMP** connector may permanently damage the amplifier and/or preamplifier. TENNELEC assumes no liability for such instrument damage.**

⊗ ⊗

TABLE OF CONTENTS

Section	Page
1.0 INTRODUCTION	1
2.0 SPECIFICATIONS	2
2.1 PERFORMANCE	2
2.1.1 SPECTROSCOPY AMPLIFIER	2
2.1.2 TIMING AMPLIFIER	3
2.2 CONTROLS	3
2.2.1 FRONT-PANEL CONTROLS	3
2.2.1.1 SPECTROSCOPY AMPLIFIER	3
2.2.1.2 TIMING AMPLIFIER	4
2.2.2 INTERNAL CONTROLS , ,	4
2.2.2.1 SPECTROSCOPY AMPLIFIER	4
2.2.2.2 TIMING AMPLIFIER	5
2.3 CONNECTORS	5
2.3.1 FRONT-PANEL CONNECTORS	5
2.3.1.1 COMMON	5
2.3.i.2 SPECTROSCOPY AMPLIFIER	5
2.3.1.3 TIMING AMPLIFIER	6
2.3.2 REAR-PANEL CONNECTORS ,	6
2.3.2.1 COMMON	6
2.3.2.2 SPECTROSCOPY AMPLIFIER	7
2.4 POWER REQUIREMENTS	7
2.5 OTHER INFORMATION	7

TABLE OF CONTENTS

Section	Page
3.0 INSTALLATION	8
3.1 POWER CONNECTION	8
3.2 PREAMPLIFIER CONNECTION	8
3.2.1 ORIGINAL TENNELEC PREAMPLIFIERS (Serial Number <2000)	8
3.2.2 TENNELEC PREAMPLIFIERS (Serial Number >2000)	9
3.2.3 OTHER PREAMPLIFIERS	9
3.2.4 GENERAL PRECAUTIONS	10
3.2.5 HIGH VOLTAGE DETECTOR BIAS	10
4.0 OPERATING PROCEDURES	10
4.1 FIRST-TIME OPERATIONS	10
4.1.1 EQUIPMENT REQUIRED	10
4.1.2 TEST SYSTEM SETUP	11
4.1.3 TEST SYSTEM OPERATION	11
4.2 POLE-ZERO (P/Z) CANCELLATION	14
4.3 AMPLIFIER NOISE	15
4.4 PEAKING TIME CONSIDERATIONS	16
4.4.1 GAS PROPORTIONAL DETECTOR	16
4.4.2 GERMANIUM DETECTOR	16
4.4.3 SCINTILLATION DETECTOR	19
4.4.4 SILICON SURFACE BARRIER DETECTOR	19
4.5 OUTPUT TERMINATION CONSIDERATIONS	19
5.0 CIRCUIT DESCRIPTION	20
6.0 STRIPPING DAMAGE	22
7.0 SERVICING	22
8.0 *****	22

TABLE OF FIGURES

Figure	Page
1.1 Typical Resolution vs Countrate	1
3.1 Original TENNELEC Preamplifier Power Connections . .	8
4.1 TC 248 Slow Linear Output Pulse Shapes	12
4.2 TC 248 Fast Timing Output Pulse Shape	13
4.3 P/E Control Adjustment	14
4.4 Resistor Bridge/Diode Limiter	15
4.5 Unipolar Noise vs Gain and Peaking Time	17
4.6 Bipolar Noise vs Gain and Peaking Time	18
4.1 Resolution vs Countrate	19
TC 248 BLOCK DIAGRAM	24
TC 248 SCHEMATIC DIAGRAM (Sheet 1 of 3)	25
TC 248 SCHEMATIC DIAGRAM (Sheet 2 of 3)	26
TC 248 SCHEMATIC DIAGRAM (Sheet 3 of 3)	27
ADT 1 PREAMP POWER ADAPTER.	28

1.0

INTROWCTION

The TENNELEC Node1 **TC 246** combines a high **perf**ormance spectroscopy **(SLOW) amplifier** with a general purpose timing **(FAST) amplifier in a single-width** NIM module. The **FAST/SLOW AMPLIFIER** is ideally suited for applications **which** require both timing and energy information from the detector: **signals.**

The **SLOW AMPLIFIER** incorporates a five-pole (two complex-conjugate-pair and one real) active **filter** network to generate a **near-gaussian** shaped unipolar pulse which **optimizes** the signal-to-noise ratio while allowing excellent count rate capability. **Switch-**selectable peaking times of 1, 3 and 6 **usec** allows the **TC 248** spectroscopy amplifier to be matched to the signal processing **requirements** of most detectors.

A **gated** baseline restorer **(BLR)** with automatic threshold and restorer rate **circuits** provides superb high count rate performance without degrading ultimate resolution.

The excellent dc stability of the **TC 248 spectroscopy** amplifier eliminates **spectrum** broadening caused by dc shift **of** the amplifier output.

A unipolar delay option **allows the** slow unipolar output to be delayed 2 **usec** for gating applications. A convenient slide-switch allows a normal (prompt) unipolar output when **the delay** is not required.

The **FAST AMPLIFIER** incorporates CR-RC shaping to **produce a high initial slope unipolar** pulse for accurate timing measurements **at low constant fractions.** Independent differentiator and **integrator** time constant switches allow the pulse decay time and **rissetime** to be custom-tailored to the experiment.,

2.0 SPECIFICATIONS

2.1 PERFORMANCE

2.1.1 SPECTROSCOPY AMPLIFIER

PULSE SHAPING Active shaping networks produce pseudo-gaussian shaped unipolar pulses with selectable peaking times* (t_p) of 1, 3 or 6 usec. Unipolar pulse width at the 1.0% level equals $2.8 t_p$. Bipolar pulse peaks* at $0.78 t_p$ and crossover* occurs at $1.36 t_p$. Bipolar crossover is delayed by $0.32 t_p$ from the prompt unipolar peak;

GAIN RANGE Continuously variable from x5 to x750.

INTEGRAL NONLINEARITY $\leq \pm 0.05\%$ over 0 to +10V output range for 3 usec peaking time.

NOISE Less than, 3.0 uV rms referred to the input for 6 usec peaking time, unipolar shaping and maximum gain; typically less than 3.2 uV rms for gain greater than 100. Bipolar typically less than 5.0 uV rms for gain greater than 100.

TEMPERATURE INSTABILITY

UNIPOLAR

Gain $\leq \pm 0.01\%/^{\circ}\text{C}$, 0 to 50°C
DC Level $\leq \pm 10 \text{ uV}/^{\circ}\text{C}$, 0 to 50°C

BIPOLAR

Gain $\leq \pm 0.01\%/^{\circ}\text{C}$, 0 to 50°C
UC Level $\leq \pm 30 \text{ uV}/^{\circ}\text{C}$, 0 to 50°C

WALK $\leq \pm 3 \text{ nsec}$ over a 50:1 dynamic range for 3 usec peaking time.

OVERLOAD RECOVERY Unipolar output recovers to within 2% of rated output from x300 overload in less than 2.5 non-overloaded pulse widths at maximum gain. Bipolar output recovers to within 2% of rated output from x300 overload in less than 2.0 non-overloaded pulse widths at maximum gain.

*Measured from the leading 1.0% of maximum signal amplitude.

SPECTRUM BROADENING (Unipolar)** Typically less than **10% FWHM** for a ^{60}Co 1.33 MeV gamma line at 85% of rated output, 3 **usec** peaking time and incoming rate of 1 to 100 kcps.

SPECTRUM SHIFT (Unipolar)** Peak position shifts typically less than **0.02%** for a ^{60}Co 1.33 MeV gamma line at 85% of rated output, 3 **usec** peaking time and incoming rate of 1 to 100 kcps.

2.1.2 TIRING AMPLIFIER

PULSE SHAPING CR-RC shaping networks produce high initial slope unipolar pulses with independently selectable **differentiation and integration** time constants of 5, 10, 20, 50 and 100 nsec.

RISETIME -10 nsec with minimum integration and maximum differentiation.

GAIN RANGE Variable from 5 to 250 in a 1-2.5-5 sequence.

TEMPERATURE INSTABILITY

Gain $\leq \pm 0.1\%/^{\circ}\text{C}$, 0 to 50°C

OUTPUT AMPLITUDE RANGE 0 to -5 volts linear into a 50 ohm load.

2.2 CONTROLS

2.2.1 FRONT-PANEL CONTROLS

2.2.1.1 SPECTROSCOPY AMPLIFIER

COARSE GAIN Six-position rotary switch selects gain factors from 10 to 500 in a 1-2-5 sequence.

FINE GAIN Ten-turn precision potentiometer with linear calibration from 500 to 1500. The dial numbers should be **considered** as **multipliers operating on the** COARSE GAIN setting, with 500 and 1500 corresponding to multipliers of 0.5 and 1.5 respectively. The FINE GAIN control extends the **total gain range** from 5 to 750.

**Results may not be reproducible if measurements are made with a detector which exhibits a larger number of slow-risetime signal components.

POLE-XERO (P/E) 15-turn **screwdriver adjustable** control for cancellation of **preamplifier decay** times from 35 **usec** to **infinity**.

POS-NEG Two-position toggle switch selects either positive or negative polarity preamplifier **signals** for the amplifier input.

BLR-P/Z **Two-position** toggle **switch** enables the **baseline** restorer in the BLR position. The P/E position **disables the baseline** restorer for accurate pole-zero cancellation adjustment.

2.2.1.2 TIWING AWPLIFIER

GAIN Six-position **rotary** switch selects gain factors from 5 to 250 in a 1-2.5-5 **sequence**.

2.2.2 INTERNAL CONTROLS

2.2.2.1 SPECTROSCOPY AWPLIFIER

PEAKING TIME **Four** individual **three-position** elide switches (accessible through the **side shield**) select unipolar peaking **time of 1, 3 or 6 usec**. Bipolar **peaking** time equals 0.78 of selected unipolar peaking **time**.

NOTE: All four PEAKING TIME switches MUST be set to the **same** position for proper operation of the amplifier.

UNIPOLAR DELAY OPTION **When installed**, a two-position slide **switch** (accessible through the side shield) selects either prompt (OUT) or a 2 **usec** delay (IN) for the unipolar output. **The bipolar** output is unaffected by the unipolar delay **selection**.

OUTPUT **IMPEDANCE** Individual **two-position** slide **switches** (accessible through the side shield) **selects an output impedance of either <1 ohm or 50 ohms** for the **UNIPOLAR** and **BIPOLAR** outputs. The **total peak** output current is limited to 50 mA, for either the **UNIPOLAR** and **BIPOLAR** outputs, when operated in the **<1 ohm** position. **All outputs** are short circuit protected in either **OUTPUT IMPEDANCE** position.

2.2.2.2 TIMING AMPLIFIER

INPUT POLARITY Two-position slide switch (accessible through the side shield) selects either **positive or negative input polarity for the timing amplifier.**

DIFFERENTIATOR TIME CONSTANT Four individual two-position slide switches (accessible through the side shield) select **differentiator time constants** of 10, 20, 50 and 100 nsec. Switches may be **used** in combination to provide intermediate values, **i.e., 50 and 100 nsec IN** to provide a 150 nsec time constant. All switches **OUT provide a 5 nsec time constant.**

INTEGRATOR TIME CONSTANT Four individual two-position slide switches (accessible **through the side shield**) select integrator time constants of 10, 20, 50 and 100 nsec. Switches **may be used** in combination to provide intermediate values, **i.e., 50 and 100 nsec IN** to provide a 150 nsec time constant. All switches **OUT provides a 5 nsec time constant.** The **rissetime** is approximately 2.2 times the selected **time constant.**

2.3 CONNECTORS

2.3.1 FRONT-PANEL CONNECTORS

2.3.1.1 COMMON

INPUT BNC connector accepts either positive or negative **polarity input signals** which are applied to both the **spectroscopy and timing amplifiers.** The input signal **rissetimes** should be **less** than the selected timing channel integrator **time constant** and the decay time should **be within the range of 35 usec to infinity.** The input **impedance is 50 ohms dc-coupled,** and the input is protected to **f3.5 volts dc absolute maximum.**

2.3.1.2 SPECTROSCOPY ANPLIFIER

UNIPOLAR BNC connector provides **,&-restored (0 ±5 mV) unipolar shaped output pulses with a full-scale range of 0 to +10V (+11.5V maximum).** The **UNIPOLAR OUTPUT IMPEDANCE** switch selects an output impedance of either **<1 ohm or 50 ohms** for both front- and rear-panel **UNIPOLAR output connectors.** The total peak output current **is limited to 50 mA** in the **<1 ohm position.** Both **UNIPOLAR** connector, **6 are short circuit protected** in either **UNIWLAR OUTPUT IWPEDANCE position.**

BIPOLAR BNC connector provides ground referenced (0 \pm 10 mV) bipolar shaped output pulses with a full-scale range of 0 to +10V (+11.5V maximum). The BIPOLAR OUTPUT IMPEDANCE switch selects an output impedance of either \leq 1 ohm or 50 ohms for both front- and rear-panel BIPOLAR output connectors. The total peak output current is limited to 50 mA in the \leq 1 ohm position. Both BIPOLAR connectors are short circuit protected in either BIPOLAR OUTPUT IMPEDANCE position.

2.3.1.3 TIMING AMPLIFIER

OUTPUT Lemo type connector provides shaped timing output pulses with a full-scale linear range of 0 to -5 volts (-7 volts maximum) when terminated into a, 50 ohm load. Risetime and decay time selected by the INTEGRATOR and DIFFERENTIATOR switch settings.


2.3.2 REAR-PANEL CONNECTORS

2.3.2.1 COMMON

INPUT Refer to INPUT of Sec. 2.3.1.1.

PREAMP A O-pin Amphenol 17-10090 type connector accepts signals from TENNELEC preamplifiers and provides power to TENNELEC or other commercially available preamplifiers.

• * • • CAUTION • • • *

The **PREAMP** connector on this instrument is directly compatible only with TENNELEC preamplifiers with serial numbers greater than 2000. It is also directly compatible with standard Aptec,  EGCG Ortec and FGT preamplifiers.

If a TENNELEC preamplifier, with serial number less than 2000 is used, then a model ADT 1 **PREAMP POWER ADAPTER** must be used (See Sec. 3.2.1).

If there are any questions regarding the compatibility of the **PREAMP** connector of this instrument, please contact the TENNELEC Marketing Department for assistance.

***** **WARNING** *****

Improper connection to the **PREAMP** connector may permanently **damage the** amplifier and/or preamplifier. **TENNELEC** assumes no liability for such instrument

PREAMPLIFIER PIN ASSIGNMENT

PIN	:	DESCRIPTION
1	:	POWER GND
2	:	SIGNAL GRD
3	:	PREAMP OUT
4	:	+12V
5	:	TEST
6	:	-24V
7	:	+24V
8	:	PHANTOM
9	:	-12V

2.3.2.2 SPECTROSCOPY AMPLIFIER

UNIPOLAR Refer to **UNIPOLAR** of Sec. 2.3.1.2.

BIPOLAR Refer to **BIPOLAR** of Sec. 2.3.1.2.

2.4 POWER REQUIREMENTS (excluding preamplifier)

+24V, 45 mA; +12V, 135 mA
 -24V, 45 mA; -12V, 115 mA

2.5 OTHER INFORMATION

OPERATING TEMPERATURE 0 to 50°C

WEIGHT Shipping 5.0 lbs (2.3 kg)
 (Net) 2.5 lbs (1.2 kg)

DIMENSIONS Standard single-wide NIM module (1.35 x 8.714 in.) per TID 20893 (Rev).

INSTRUCTION MANUAL One provided with each instrument ordered.

WARRANTY One year.

3.9 INSTALLATION

3.1 POWER CONNECTION

The TC 248 FAST/SLOW, AMPLIFIER requires a NIM-standard bin and power supply, such as the TENNELEC TB3/TC 911, for operation. The bin provides mechanical mounting and power supply distribution. Always turn OFF the bin power supply when inserting or removing any modules.

All TENNELEC NIM modules are designed so that it is not possible to overload the power supply, even with a full complement of modules in the bin. Since this may not be true when the bin contains modules other than those of TENNELEC design, the power supply voltages should be checked after all modules have been inserted. The TENNELEC Bin and Power Supply provides power supply test points on the bin control panel for monitoring the dc voltage levels.

3.2 PREAMPLIFIER CONNECTION

3.2.1 ORIGINAL TENNELEC PREAMPLIFIERS (Serial Number <2000)

The **PREAMP** connector of this amplifier is NOT directly compatible with TENNELEC preamplifiers with serial numbers less than 2000. An optional model ADT 1 PREAMP POWER ADAPTER is required.




Fig. 3.1 Original TENNELEC Preamplifier Power Connections

The ADT 1 PREAUP POWER ADAPTER should be physically mounted to the amplifier's rear-panel **PREAMP** connector (See Fig. 3.1).

The preamplifier's power cable is then connected to the ADT 1 adapter. As with the original TENNELEC configuration, the preamplifier's output signal is present in the power cable and no external BNC cable is required from the preamplifier's output to the amplifier's input. The amplifier's front-panel PDS-NEG switch should be set to match the preamplifier's output signal polarity.

3.2.2 TRNNELBC PREAMPLIFIERS (Serial Number > 2000)

The **PREAMP** connector of this amplifier is directly compatible with TENNELEC preamplifiers with serial numbers greater than 2000.

The preamplifier's power cable is directly connected to the amplifier's **PREAMP** connector. The preamplifier's output signal is present in the power cable and no external BNC cable is required from the preamplifier's output to the amplifier's input. The amplifier's front-panel PGS-NEG switch should be set to match the preamplifier's output signal polarity.

3.2.3 OTHER PRBAWPLIFIERS

The **PREAMP** connector of this amplifier is directly compatible with standard Aptec, Canberra, EG&G Ortec and PGT preamplifiers.

The preamplifier's power cable is directly connected to the amplifier's **PREAMP** connector. The preamplifier's output signal should be connected to the amplifier's rear-panel INPUT connector with a BNC cable. The BNC cable should be spirally wrapped with the preamp power cable to reduce noise pickup.

3.2.4 GENERAL PRECAUTIONS

When an external BNC cable longer than ten feet is used to connect the preamplifier output to the amplifier input, the characteristic impedance of the cable should match the impedance of the preamplifier output. All TENNELEC preamplifiers contain 50 ohm series termination, therefore 50 ohm RG-58 cable is recommended.

To minimize electrical noise pickup when external BNC signal and/or high voltage cables are connected between the preamplifier and NIM bin (containing the amplifier and high voltage power supply), the cable lengths should match the preamp power cable length and be spirally wrapped with the preamp power cable.

3.2.5 RICE VOLTAGE DETECTOR BIAS

It is recommended that the detector high-voltage power supply be mounted in the same NIM bin as the amplifier to reduce ground-loop noise pickup. The high-voltage cable between the H.V. supply and the preamplifier should be spirally wrapped with the preamp power cable.

4.0 OPERATING PROCEDURES

4.1 FIRST-TIME OPERATIONS

Users will find it helpful to familiarize themselves with the TC 248 FAST/SLOW AMPLIFIER by conducting a few simple tests.

4.1.1 EQUIPMENT REQUIRED

1. NIM Bin and Power Supply (TENNELEC TB3/TC 911 or equivalent).
2. Precision Tail Pulser (TENNELEC TC 812 or equivalent).
3. Oscilloscope (TEKTRONIX 465 or equivalent).
4. Digital Multimeter (FLUKE 8010A or equivalent).
5. Shielded 50 ohm cables with BNC connectors.
6. BNC tee and 50 ohm terminator.
7. BNC to Lemo Adapter.

4.1.2 TEST SYSTEM SETUP

Set the TC 240 controls as follows, then insert into the NIM bin such that the internal controls are accessible.

COARSE GAIN	20
FINE GAIN	1.00
P/Z	Fully CCW
POS-NEG	POS
BLR-P/Z	P/Z
TIMING AMP GAIN	10
PEAKING TIME (PCB)	3 usec
UNIPOLAR DELAY (PCB)	OUT
TIMING AMP INPUT POLARITY (PCB)	POS
TIMING AMP DIFFERENTIATOR (PCB)	100 nsec
TIMING AMP INTEGRATOR (PCB)	100 nsec

Set the Tail Pulser controls as follows:

PULSE HEIGHT	10.0
RELAY	OFF
POL	+
DIRECT OUT/EXT IN	DIRECT OUT
OUTPUT	5V
ATTENUATION	10

Set the Oscilloscope controls as follows:

CH A VERT SENS	5 Volts/Div (dc-coupled)
CH B VERT SENS	5 Volts/Div (dc-coupled)
VERT DISPLAY MODE	CH A
HORIZ SWEEP	2 usec/Div
TRIGGERING	POS EXTERNAL

Connect the DIRECT OUTPUT of the TC 812 to the EXTERNAL TRIGGER INPUT of the oscilloscope with a BNC cable. Connect a BNC cable between the TC 812 OUTPUT and the INPUT of the TC 248 amplifier.

4.1.3 TEST SYSTEM OPERATION

Apply power to the NIM bin and set the Digital Multimeter (DMM) for 200 mV dc full-scale. Connect the DMM's negative lead to the GND (ground) test point of the NIM bin. Touch the DMM's positive lead to the center pin of the TC 248's UNIPOLAR connector; the DMM should indicate 0 \pm 5 mV. Repeat for BIPOLAR connector; the DMM should indicate 0 \pm 10 mV.

Connect BNC cables from the UNIPOLAR and BIPOLAR outputs of the TC 248 to the Channel A and Channel B vertical inputs, respectively, of the oscilloscope. Turn the RELAY power of the TC 812 to the ON position and adjust the oscilloscope trigger control for a stable display. Adjust the TC 248 P/Z control for proper compensation (refer to Sec. 4.2). Return the oscilloscope vertical and horizontal controls to their initial settings. Switch the oscilloscope to the alternate display mode. Positive unipolar and bipolar shaped pulses should be present (See Fig. 4.1). Note the time relationship of the unipolar peak to both the bipolar peak and bipolar zero-crossing. Switch the TC 248 BLR-P/Z control to the BLR position; the unipolar signal should not change.




Fig. 4.1 TC 248 Slow Linear Output Pulse Shapes

For instruments with the UNIPOLAR DELAY OPTION, position the unipolar signal peak around the center vertical graticule line. Switch the UNIPOLAR DELAY control to the IN position. Note that the unipolar signal is delayed by 2 usec, while the bipolar signal remains prompt. Return the UNIPOLAR DELAY control to the OUT position.

Change all four of the TC 248 PEAKING TIME switches to the 6 usec position. The unipolar signal should now peak at 6 usec.

Reduce the TC248's **COARSE GAIN** to 10. The maximum amplitude of both the **unipolar** and **bipolar** signals should decrease to 5.0 volts. Increase the TC 248's **FINE GAIN** control to 1.500 (fully CW). The maximum amplitude of both signals should increase to 7.5 volts.

Remove both **cables** between the **amplifier's** outputs and the **oscilloscope's** inputs. Place a **BNC tee** on the **CH B** input of the **oscilloscope** and connect a 50 ohm terminator on one side of the tee. Connect a **BNC cable** to the remaining side of the **BNC tee**. Attach a **BNC-to-Lemo adapter** to the remaining end of the **BNC cable** and connect to the TC 248 **TIMING AMP OUTPUT**.

Set the Oscilloscope controls as follows:

CR B VERT SENS	0.5 Volt/Div (dc-coupled)
VRRT DISPLAY MODE	CH B
HORIZ SWEEP	0.2 usec/Div
TRIGGERING	POS EXTERNAL

Negative 2.0 volt **fast-shaped** pulses should be present (See Fig. 4.2). Reduce the TC 248's **TIRING AWP GAIN** to 5. The maximum amplitude of the timing signal should decrease to 1.0 volts.




Fig. 4.2 TC 248 Fast Timing Output Pulee Shape

4.2 POLE-ZERO (P/Z) CANCELLATION

Accurate setting of the P/Z control is essential for good resolution at high count rates in unipolar operation and for correct operation of the BLR circuit. With bipolar operation, accurate setting is not important regarding resolution, but is important if quick recovery from heavily overloading signals is required. The adjustment procedure is as follows:

Using a detector, and a radioactive source as an input, observe the UNIPOLAR output signals of the amplifier on an oscilloscope with a triggered sweep. The oscilloscope **MUST** be dc-coupled to the amplifier.




Fig. 4.3 P/Z Control Adjustment

Adjust the amplifier gain so that the highest energy peak produces an amplifier output of approximately 9 volts. Adjust the source strength and/or spacing from detector to provide a count rate between 2,000 and 10,000 cpe. Increase the oscilloscope vertical sensitivity to 100 mV/Div.

Compensation is accomplished by turning the P/2 control for flattest possible baseline (See Fig. 4.3).

A second method of setting the control is to increase the sweep speed 80 that the pulse duration occupies approximately 2 cm of sweep. Then; looking at the baeeline (which will be fuzzy) about one pulse width from the end of the pulse, adjust the P/K control for minimum baseline smear.

NOTE: Oscilloscopes such as the TBKTRONIX Model 465 and 475 will overload with a 10 volt input signal when the vertical sensitivity is eet for 100 mV/Div or less. The Resistor Bridge/Diode Limiter shown in Fig. 4.4 is recommended to prevent overloading the oscilloscope (this bridge/limiter is also used for amplifier linearity measurement),,, The FROM GENERATOR OUTPUT connector is not uied when connected as an oscilloscope overload limiter.




Fig. 4.4 Resistor Bridge/Diode Limiter

4.3

ARPLIFIBRNOISE

The typical equivalent noise referred to the input vs amplifier gain for unipolar and: bipolar rhapsing is shown in Fig. 4.5 and Fig. 4.6 respectively.




Fig. 4.5 UNIPOLAR NOISE vs GAIN and PEAKING TIME




Fig. 4.6 BIPOLAR NOISE vs GAIN and PEAKING TIME

Lowest equivalent input noise results when a high COARSE GAIN and low FINE GAIN are used to produce the same overall gain. For example, with unipolar shaping and a 3 usec peaking time, a COARSE GAIN of 100 and a FINE GAIN of 0.750 result in an overall gain of 75 with an equivalent input noise of 4.8 uV, while a COARSE GAIN of 50 and a FINE GAIN of 1.500, although still producing an overall gain of 75, result in an equivalent input noise of 8.2 uV.

4.4 PEAKING TIME CONSIDERATIONS

The optimum peaking time for a particular system depends on the detector characteristics and counting rate. A general discussion of peaking time requirements is presented below; consult the detector manufacturer for specific shaping requirements.

4.4.1 GAS PROPORTIONAL DETECTOR

The required peaking time for a gas proportional detector depends on the charge collection time of the detector, which is related to the physical size (both the outer electrode and center wire), fill gas and high voltage. Generally, a larger detector will have a slower collection time.

Correct pole-zero compensation is not possible because of the peculiar charge collection characteristics, resulting in either an undershoot or back porch on the trailing-edge of a unipolar shaped signal. Bipolar shaping is recommended because it reduces this effect without degrading detector resolution.

4.4.2 GERMANIUM DETECTOR

Unipolar shaping and a peaking time of 3 or 6 usec generally will give best results for germanium detectors, depending on the detector size, configuration, and counting rate.

A peaking time of 6 usec is preferred for applications requiring ultimate low-count rate resolution, while 3 usec is preferred at high-count rates. Resolution as a function of count rate at 3 and 6 usec peaking times for a typical intrinsic germanium detector is shown in Fig. 4.7. A performance crossover occurs above 80 Kcps with 3 usec peaking time providing better resolution.




Fig. 4.7 Resolution vs Count Rate

4.4.3 SCINTILLATION DETECTOR

Scintillation detectors require a peaking time of two to three times the decay time constant of the scintillator (a 1 usec peaking time is about optimum for NaI scintillators). Bipolar shaping is preferred because it reduces the effect of overload and detector microphonics without degrading resolution.

4.4.4 SILICON SURFACE BARRIER DETECTOR

Unipolar shaping and a peaking time of 1 or 3 usec usually will give best results, depending on the detector capacitance and leakage current. Generally, large area detectors are resolution-limited by their high capacitance and thus benefit from a longer peaking time. Small detectors are resolution-limited by their leakage current, which implies a shorter peaking time.

4.5 OUTPUT TERMINATION CONSIDERATIONS

The TC 248 allows individually switchable output terminations of either <1 ohm or 50 ohms for both the UNIPOLAR and BIPOLAR outputs. The TC 248 is shipped with both OUTPUT IMPEDANCE switches in the 50 ohm position.

Series 50-ohm, termination is recommended for general applications and is required for long cable lengths and/or high noise environments. In extreme situations series termination along with shunt termination at the receiving end may be required, however, this reduces the signal amplitude at the receiving end of the cable to .501 of the non-shunt-termination value. An in-line 50, ohm terminator or BNC Tee and 50 ohm terminator may be used for shunt termination.

The <1 ohm termination is useful when driving several instruments from the same amplifier "output connector. The output signal amplitude will be relatively independent of the load impedance, which is not the case with series 50 ohm termination. However, the amplifier may oscillate when driven into overload.

5.0

CIRCUIT DESCRIPTION

The TC 248 FAST/SLOW AMPLIFIER module contains a high performance spectroscopy amplifier and a general purpose timing amplifier which share a common input.

The spectroscopy amplifier consist of three gain stages, two active-filter stages, two output driver stages, and an active, baseline restorer, (BLR) with automatic threshold and restorer rate, circuits as shown in the. TC 248 BLOCK DIAGRAM. The overall transfer function of the amplifier synthesize an $e^{-3t} \sin^4 t$ time response resulting in a pseudo-gaussian pulse shape.

Transistors Q1 through Q6 form a low-noise input op-amp stage with switchselectable gains of 20, 10, 4, or 2. Diodes D1 and D2 provide input protection, while diodes D3 through D5 provide biasing and overload limiting.

Input-polarity selection and additional gain is provided by IC1. Diode D8 provides overload limiting.

Fine-gain adjustment and additional switch selectable gain is provided by IC2. Diode D9 Provides overload limiting.

IC3 and IC1 form a four-pole active integrator.:

The signal from the integrators is differentiated, then buffered by IC5 and transistors Q8 and Q9 to provide a high-current BIPOLAR output. Diodes D13 and D14 provide overload limiting.

The integrated signal is directly applied to IC6 and transistors Q10 and Q11 to provide a high-current UNIPOLAR output. Diodes D20 and D21 provide overload limiting.

Baseline restoration of the UNIPOLAR output is provided by transconductance amplifier IC7. Capacitor C51 functions as an output level memory. IC9 and IC10 form a peak detector for the negative noise excursions at the UNIPOLAR output, which determines the BLR gate threshold. IC8 and transistor Q13 gate the restorer current whenever the UNIPOLAR output exceeds the BLR gate threshold. An additional, exponentially decaying restoration current is provided when switch S9 is in the BLR position. This additional current is disabled with switch S9 in the P/Z position for accurate pole-zero cancellation adjustment.

The timing amplifier consists of an input buffer, monolithic amplifier with differential outputs, and gain of x10 output driver stage.

Transistors Q15 through Q18 form a high input impedance, low output impedance buffer amplifier which isolates the timing amplifier differentiator from the input signal.

Reed Relays RLY1 through RLY4 select an input attenuation factor of x1, x2, x5 or x10 depending on the setting of the timing amplifier GAIN switch.

Monolithic video amplifier - IC11 amplifies the differentiated signal and provides both inverted and non-inverted outputs. Relays RLY5 and RLY6 select the gain of the video amplifier IC by changing the emitter resistance of the differential input stage.

Input polarity selection is accomplished by choosing the proper output signal from the video amplifier.

Transistors Q19 through Q22 form a gain of x10 non-inverting output driver stage.

6.0 SHIPPING DAMAGE

Upon receipt of the instrument, examine it for shipping damage. Damage claims "should be filed with the carrier. The claims agent should receive a full report: a copy of that report should be sent to TENNELEC, Inc., P.O. Box 2560, Oak Ridge, Tennessee 37831-2560. The model number and serial number of the instrument must be included in the report. Any remedial action taken by TENNELEC, Inc., will be based on the information contained in this report.

7.0 SERVICING


In the event of a component failure, replacement may be done in the field or the instrument may be returned to our plant for repair. There will be no charge for repairs that fall within the warranty period.

8.0 WARRANTY

In connection with TENNELEC's warranty (inside front cover), TENNELEC suggests that if a fault develops, the customer should immediately notify the TENNELEC Customer Service Manager. He may be able to prescribe repairs and send replacement parts which will enable you to get the instrument operating sooner and at less expense than if you returned it.

Should return prove necessary, the TENNELEC Customer Service Manager must be informed in WRITING, BY CABLE or TWX of the nature of the fault and the model number and serial number of the instrument. Pack the instrument well and ship PREPAID and INSURED to TENNELEC, Inc., 601 Oak Ridge Turnpike, Oak Ridge, Tennessee 37830. As stated in the warranty DAMAGE IN TRANSIT WILL BE REPAIRED AT THE SENDER'S EXPENSE as will damage that obviously resulted from abuse or misuse of the instrument.



Quotations for repair of such damage will be sent for your approval before repair is undertaken.



NOTES:
 (1) * DENOTES SWITCHED 'PEAKING TIME' COMPONENT.
 ** DENOTES SWITCHED FAST AMP 'SHAPING CONSTANT' COMPONENT.
 (2) O DENOTES TEST POINTS.
 (3) TYPICAL WAVEFORMS SHOWN ARE FOR THE FOLLOWING CONDITIONS:
 (A) COARSE GAIN: 10
 (B) FINE GAIN: 1.000
 (C) POS-NEG: NEG
 (D) BLR - P/2: P/2
 (E) TIMING AMP GAIN: 5
 (F) -1VOLT STEP INPUT

REV	ECN	DESCRIPTION	BY	DATE
0		RELEASE TO PRODUCTION	DA/SMH	8-27-64
REVISIONS				
TENNELEC OAK RIDGE, TENNESSEE				
TC 248 BLOCK DIAGRAM				
<small> DESIGNED BY: D. FOSSELL, I-2034 DRAWN BY: S. HINSHAW CHECKED BY: J. W. ... APPROVED BY: ... PRODUCT NO: 1400-1187 </small>				
				TC 248

REVISIONS					
REV	ECN	DESCRIPTION	INC BY	APP BY	DATE
0		RELEASED TO PRODUCTION	D.E.	SM	6-27-59




TENNELEC
 OAK RIDGE, TENNESSEE

TC 248
 SCHEMATIC DIAGRAM

DESIGNED BY	DATE	REV	BY	CHK
D. GOSSEL	1-25-59	1	S. MINSHAW	SM
APPROVED BY				
A. S. GIBSON	5-24-59		S. MINSHAW	SM
ISSUED BY				
				TC 248

100-1187

REVISIONS					
REV	ECN	DESCRIPTION	HC BY	APP BY	DATE
0		RELEASED TO PRODUCTION	D.E.	SMAN	6-27-84



- NOTES:
- (1) ALL RESISTORS IN OHMS UNLESS OTHERWISE NOTED.
 - (2) ALL RESISTORS 1/4 WATT 5% D.C. UNLESS OTHERWISE NOTED:
 - 1/4 WATT 1% METAL FILM (RN 55C).
 - 1/8 WATT 5% D.C.
 - 1/8 WATT 1% METAL FILM (RN 50C).
 - 1/4 WATT 5% CARBON COMP.
 - (3) ALL CAPACITORS ARE IN MICROFARADS UNLESS OTHERWISE NOTED.
 - (4) ALL CAPACITORS ARE CERAMIC UNLESS OTHERWISE NOTED:
 - 1% DIPPED MICA.
 - TANTALUM.
 - MYLAR.
 - (5) O DENOTES TEST POINTS.
 - (6) ALL NPN TRANSISTORS ARE 2N3904 UNLESS OTHERWISE NOTED.
 - (7) ALL PNP TRANSISTORS ARE 2N3906 UNLESS OTHERWISE NOTED.
 - (8) ALL DIODES ARE 1N4154 UNLESS OTHERWISE NOTED.

LAST USED	NOT USED
R190	
C107	
D52	
Q22	
IC11	
S14	
L10	
CN9	
TR20	
DL1	


TENNELEC
OAK RIDGE, TENNESSEE

TC 248
SCHEMATIC DIAGRAM

DESIGNED BY	D. E. SELLS	DATE	7-31-84	REVISED BY	S. HINSHAW	DATE	8-1-84	REV. NO.	3	OF	3
APPROVED BY	A. J. ...	DATE		APPROVED BY		DATE		PROJECT NO.	1400-1187		

TC 248

(1) ADT-1 ADAPTER PHYSICALLY MOUNTED
AT AMPLIFIER **PREAMP** CONNECTOR.



PREAMPLIFIER PIN ASSIGNMENT

PIN	ORIGINAL TENNELEC	NEW TENNELEC
1	PREAMP OUT	POWER GND
2	PHANTOM	SIGNAL GND
3	+12V	PREAMP OUT
4	TEST	+12V
5	HV BIAS	TEST
6	SIGNAL GND	-24V
7	POWER GND	+24V
8	-12V	PHANTOM
9	TEST GND	-12V

QTY	ITEM	PART NO.	DESCRIPTION	REMARKS
LIST OF MATERIALS OR PARTS LIST				
UNLESS OTHERWISE SPECIFIED, DIMENSIONS ARE IN INCHES		TENNELEC		
TOLERANCES: DECIMALS .XXX ± .005 .XX ± .015 .X ± .020 FRACTIONS ANGLES ±1/64 ±1°		P. O. BOX D, OAK RIDGE, TENNESSEE 37830		
		PREAMPLIFIER POWER ADAPTER		
SIZE	MATERIAL	DRN BY	DATE	ENGINEER
⌀		D. EDSSELL	1-7-81	S. HINSHAW
		DATE	THRO APP	DATE
		W. C. BYRNE	1-7-81	S. W. WOODS
		SCALE	DATE	DATE
		NONE		
				ADT-1