

# TENNELEC


## INSTRUCTION MANUAL

TC 307

# **INSTRUCTION MANUAL**

**TC 307 LINEAR GATE**

TC 307  
TABLE OF CONTENTS

<u>SECTION</u>	<u>PAGE</u>
1.0 INTRODUCTION . . . . .	1
2.0 SPECIFICATIONS . . . . .	1
2.1 OTHER <b>INFORMATION</b> . . . . .	4
3.0 CONTROLS. . . . .	4
3.1 FRONT-PANELCONTROLS . . . . .	.
3.2 CIRCUIT BOARD CONTROLS . . . . .	5
3.3 REAR PANEL CONTROLS . . . . .	5
4.0 CONNECTORS . . . . .	5
4.1 FRONT PANEL CONNECTORS . . . . .	5
4.2 REAR PANEL CONNECTORS . . . . .	5
5.0 FIRST TIME OPERATION . . . . .	6
5.1 <b>INSTRUMENTS</b> REQUIRED . . . . .	6
5.2 CONTROL SETTINGS . . . . .	7
5.3 WAVEFORM OBSERVATIONS, LOW-LEVEL SIGNALS . . . . .	<b>8</b>
5.4 WAVEFORM OBSERVATION, HIGH-LEVEL SIGNALS . . . . .	9
5.5 WAVEFORM OBSERVATIONS, SNITCHING SPIKES . . . . .	<b>10</b>
6.0 ADJUSTMENTS . . . . .	1 1
6.1 PEDESTAL CONTROL . . . . .	11
6.1.1 PEDESTAL ADJUSTMENT WITH SIGNAL-SOURCE OFFSET CONTROL . . . . .	11
6.1.2 PEDESTAL ADJUSTMENT, SIGNAL-SOURCE OFFSET NOT ADJUSTABLE. . . . .	<b>12</b>
6.2' INTERNAL ADJUSTMENTS . . . . .	12
7.0 CIRCUIT <b>OPERATION</b> . . . . .	15
7.1 INPUT STAGE  . . . . .	15
7.2 GATE SWITCH . . . . .	16

INTRODUCTION

The TENNELEC TC 307 Linear Gate is a single-width NIM designed to block or pass a linear signal at the command of an externally-applied gating pulse. The gain is unity, noninverting.

The TC 307 is dc-coupled, allowing it to operate from dc to high count rates without baseline shift. In the normally open (N.O.) mode, it will pass **0-10V** signals of either polarity (dc, unipolar, bipolar, sinewave, etc.) **unless** inhibited by an external gate pulse ( **$\geq +1.3V$** ). In the normally closed (N.C.) mode, it will block signals of either polarity unless enabled by the external gate pulse ( **$\geq +1.3V$** ). Gating time is  **$<15$**  ns. The gate trigger circuit remains in its **switched** state for the duration of the externally-applied gate signal and recovers within 15 ns after removal of that signal. By applying a train of gate pulses, a linear dc signal may be chopped synchronously with the gate pulses.

An LED indicates the presence of output signals (pulse or dc, either polarity) larger than **0.25V**.

Front-panel connectors are duplicated on the rear,

## TABLE OF CONTENTS

<u>SECTION</u>	<u>PAGE</u>
7.3 OUTPUT STAGE . . . . .	16
7.4 COUNTING MONITOR . . . . .	16
8.0 OPERATINGNOTES . . . . .	17
8.1 GATE CONTROL . . . . .	17
8.2 OUTPUT IMPEDANCE . . . . .	17
8.3 <1 OHM OPERATION . . . . .	17
8.4 50 OHM OPERATION . . . . .	18
8.5 PLACEMENT IN A SYSTEM . . . . .	18
8.6 AC VS DC COUPLING . . . . .	19
8.7 INCREASING <b>THE RANGE</b> OF <b>THE</b> PEDESTAL CONTROL . .	20
9.0 SHIPPING DAMAGE . . . . .	20
10.0 SERVICING . . . . .	20
<b>11.0 WARRANTY . . . . .</b>	<b>20</b>

## LIST OF FIGURES

<u>SECTION</u>	<u>PAGE</u>
<b>5-1</b> Setup for observing operation of the TC 307 . . . . .	7
<b>5-2</b> TC 307 and gate-pulse waveforms, low-level signal . .	<b>8</b>
<b>5-3</b> TC 307 and gate-pulse waveforms, high-level signal . .	9
<b>5-5</b> TC 307 switching spikes . . . . .	10
6-2 Test setup for <b>adjustment of TC 307</b> internal controls . . . . .	13
6-3 Waveforms for TC 307 setup . . . . .	14
6-4 <b>THERM BAL</b> trimmer adjustment . . . . .	15

2.0

SPECIFICATIONS

GAIN	0.99 ± 0.02
INPUT	
Coupling	DC
Polarity	Either
Range	
Linear	>10V
Maximum	
DC	±24V
Pulse	±75V, 10% duty cycle
Impedance	
Below 11V	1950 ohms ± 1%
Above 11V	460 ohms
OUTPUT,	
Coupling	DC
Polarity	Same as input (noninverting)
Range	
Open Circuit	
Linear	>10V
Maximum	±12V
Terminated, in 50 ohms	
Linear	>5V
Maximum	±6V with output R set for 50 ohms ±12V with output R set for <1ohm
Impedance	
50 ohm setting	50 ohms ±1%
<1 ohm setting	<1 ohm
PEDESTAL (Input terminated in 50 ohms)	<5mV decaying to zero in 15us
GATE SWITCHING TRANSIENTS	
Amplitude	<50 mV
Width	<20 ns full width at half max.

GATE INPUT	
Coupling	DC
Polarity	Positive
Threshold	$\leq +1.3V, 1.2V$ typical
Maximum Input	
DC	$\pm 20V$
Pulse	$\pm 75V, 10\%$ duty cycle
Hysteresis	$0.4V$ typical
Width	Gate signal is initiated when input exceeds the threshold level and remains on until input drops $0.4V$ below threshold.
Input Impedance	
Below 4V	1,000 ohms nominal
Above 4V	1,300 ohms nominal
PROPAGATION DELAY	
Linear channel	30 ns typical
Gate ( $0.1V$ overdrive)	15 ns typical
FEEDTHROUGH	
(Gate closed, $\pm 10V$ quasi-gaussian pulse, 2 us peaking time)	$\leq 5$ mV
INTEGRAL NONLINEARITY	
0 to $\pm 10V, 1k$ load, quasi-gaussian input pulse, 2us peaking time	$< 0.05\%$ referred to 10V full scale
<b>RISETIME</b> (10% - 90%)	$\leq 40$ ns, 30 ns typical
0 to $\pm 10V, 1k$ load	
<b>OVERSHOOT</b>	$< 2\%$
NOISE	
Wideband, 50 ohm input (measured with Hewlett-Packard 400 Series ac voltmeter)	200uV rms typical
Output filtered with quasi-gaussian shaping network, 2us peaking time	35uV typical

TEMPERATURE COEFFICIENT

Gain  
Zero

$\leq 100 \text{ ppm}/^{\circ}\text{C}, 0^{\circ} \text{ to } 50^{\circ}\text{C}$   
 $\leq 100 \text{ uV}/^{\circ}\text{C}, 0^{\circ} \text{ to } 50^{\circ}\text{C}$

SUPPLY VOLTAGE SENSITIVITY

Gain change  
Zero shift  
**+24V**, -24V supplies  
**+12V**, -12V supplies

$\leq 0.01\%/ \text{Volt}$   
 $\leq 0.5 \text{ mV}/ \text{Volt}$   
 $16 \text{ mV}/ \text{Volt}$

SIGNAL INDICATOR LED

Threshold

$\pm 0.25\text{V}$  typical

2.1

OTHER INFORMATION

POWER REQUIREMENTS

**+24V**  
**+12V**  
-12v  
-24V

45 mA maximum  
55 mA maximum  
50 mA maximum  
35 mA maximum

WEIGHT

Shipping  
Net

4.0 lbs. (**1.82** kg)  
2.1 lbs. (0.96 kg)

DIMENSIONS

Standard single-width NIM  
(1.35 " x 8.714") per TID  
20893 (Rev).

WARRANTY

One year

INSTRUCTION MANUAL

One provided with each  
instrument ordered

3.0

CONTROLS

3.1

FRONT-PANEL CONTROLS

PEDESTAL

Multiturn, screwdriver  
adjustable

Range

0 to  $\pm 0.25\text{V}$  nominal

MODE

Z-position toggle switch  
N.O. (normally open) and  
N.C. (normally closed).



3.2 CIRCUIT BOARD CONTROLS

N.O. ZERO	One-turn trimmer, <b>screw-</b> driver adjustable
<b>N.C.</b> ZERO	One-turn trimmer <b>screw-</b> driver adjustable
OUTPUT ZERO*	One-turn trimmer screw- driver adjustable
TBERM BAL	One-turn <b>trimmer screw-</b> driver adjustable
OUTPUT IMPEDANCE	<b>2-position</b> slide switch, 50 ohms and <b>&lt;1 ohm</b>

3.3 REAR PANEL CONTROLS None

4.0 CONNECTORS

4.1 FRONT PANEL CONNECTORS

GATE	BNC (UG 1094/U)
INPUT	BNC (UG <b>1094/U</b> )

4.2 REAR PANEL CONNECTORS

GATE	<b>BNC (UG1094/U). Parallels</b> front-panel GATE connector
INPUT	<b>BNC (UG1094/U). Parallels</b> front-panel INPUT connector
OUTPUT	<b>BNC (UG1094/U) Parallels</b> front panel OUTPUT connector when OUTPUT IMPEDANCE switch is in the <1 ohm position. Front and rear panel OUTPUT connectors are isolated <b>from</b> each other by independent 50 ohm <b>±1%</b> resistors when the OUTPUT IMPEDANCE switch is in the 50 ohm position.

---

\*Serial numbers above 209.

5.0 FIRST TIME OPERATION

The following test **is** a **good** way for users to familiarize themselves with the operation of the TC 307. The test results given below are typical.

5.1 Instruments Required

TC 307 LINEAR GATE

RECTANGULAR PULSE GENERATOR: Adjustable amplitude, **risetime 0.25us** or less, output impedance 100 ohms or less.

ATTENUATOR BOX: Impedance to match pulse generator.

OSCILLOSCOPE: 5 MHz or greater bandwidth, delayed input, externally triggered sweep, dual trace.

SHAPING AMPLIFIER: Nuclear pulse amplifier with peaking time 0.25 us **or greater.**

TERMINATOR for the pulse generator.

COAXIAL CABLES: As needed.

NIM BIN or other power source.

Connect the system as shown in Fig. 5-1 and turn it on. Note the following:

1. The TC 307 was shipped with internal adjustments accurately set for power supply voltages within **0.1V** of the nominal values of **±24V, ±12V.**
2. The OUTPUT IMPEDANCE switch was set to the **50** ohm position prior to shipment.
3. The TC 307 will not be damaged by inserting or removing it from the bin with the power ON.




Fig. 5-1 Setup for observing operation of the TC 307

## 5.2 Control Settings

With the test instruments listed below, the following control settings should be made. Settings may require modification if instruments other than those listed below are used.

OSCILLOSCOPE	TEKTRONIX 475 (200 MHz B.W.)
Channel A	0.5V/Div, dc coupled
Channel B	0.05V/Div, dc coupled
Sweep	0.5 us/Div
Tc 307	
MODE	N.C.
SHAPING AMPLIFIER	TBNNELEC TC 205A
Shaping Constant	0.25 us, Unipolar
Gain	
Coarse	10
Fine	X0.6
Polarity	+In, +Out
PULSE GENERATOR	TEKTRONIX 114
Period	0.5 ms
Width	2.7 us
Amplitude	0.9V or 1.5V. See text.
<b>ATTENUATOR</b>	X100

Waveform observations, Low-Level Signals




Fig. 5-2. TC 307 and gate-pulse waveforms, low-level signal.

In Fig. 5-2, a multiple exposure oscillogram is displayed of the gate-pulse generator output and TC 307 output. The gate was in the N.C. mode.

With the pulse generator output at **1.5V** (top waveform), the gate is opened and transmits the **90 mV quasi-gaussian** pulse from the TC 205A amplifier. Note the switching spike approximately **30 mV** high at the instant of gate closure (at **2.7 us**). Another **smaller** and narrower spike exists at the start of the pulse, but is not large enough to be visible in the figure. The delay in the onset of the TC **205A** output pulse is inherent in the TC 205A. The pedestal control was tweaked to align the baseline before and after gate closure (at **2.7 us**).

With the pulse generator output reduced to **0.9V**, the gate remains closed and no signal is transmitted. This is illustrated by the lowest line in the figure, which was displaced (with the oscilloscope positioning control) to distinguish it from the waveform above it. **The L.E.D.** should not glow under the conditions of the preceding test.

5.4 Waveform Observation, High-Level Signals

Fig. 5-3 shows the feedthrough of a **10V** signal. To obtain the waveforms, reset the controls as follows (control settings other than those listed remain unchanged):

OSCILLOSCOPE Channel B	<b>5V Div</b>
TC <b>307</b> MODE	N.O.
<b>ATTENUATOR</b>	out
PULSE GENERATOR	<b>1.0V</b>
SHAPING AMPLIFIER Fine Gain	Adjust for 10.V output <b>signa</b>




Fig. 5-3 TC 307 and gate-pulse waveforms, high-level signal.

In the figure, the center waveform (amplitude scale at the right) is that transmitted by the TC 307. The uppermost **wavetorm** is that of the pulse generator.

If the TC 307 MODE switch is set to N.C. and the oscilloscope sensitivity increased to 10 **mV/Div**, the lowest waveform **results**. (The oscilloscope position control was readjusted to separate the two lowest waveforms). This shows the feedthrough of the **10V** pulse when the gate is closed. Note that the peak of **the** feedthrough signal is synchronized with that portion of the input signal which exhibits the greatest slope.

The L.E.D. should glow when the TC 307 is transmitting the **10V** signal, but **not** when the signal transmission is blocked (gate **N.C.**).

The waveform to the right of the 3 us mark is a response of the amplifier to the trailing-edge **pulse-generator** transition and has no bearing on this test.

5.5 Waveform Observations, Switching Spikes

Disconnect the amplifier from the TC 307 and terminate its input with 50 ohms.

Change control settings as follows:

OSCILLOSCOPE	
Channel B	20 mV/Div
Sweep	0.05 us/Div

PULSE GENERATOR	
Width	0.2 us
Amplitude	<b>1.5V</b>

ATTENUATOR	X100
------------	------

The upper waveform in Fig. 5-5 is that of the pulse generator and the lower is that of the switching spikes.




Fig. 5-5. TC 307 switching spikes.

The delay between the gate-pulse transitions and the spikes is due mainly to the propagation delay ( $\approx 30$  ns) of the linear circuits in the TC 307.

In the TC 307 used for this **test**, the pedestal was about 2 mV high and caused a vertical shift of the baseline between the leading-edge and trailing-edge spikes. This pedestal was tuned out with the PEDESTAL control prior to making the oscillogram.

## 6.0 ADJUSTMENTS

### 6.1 PEDESTAL Control

In normalpoff operation, when a gate enabling pulse occurs, a transient **baseline** offset known as a "pedestal" may appear at the OUTPUT of the TC 307. Linear signals will ride on this pedestal, causing an error in signal height equal to the pedestal height. The pedestal may be adjusted to zero using the PEDESTAL control.

The PEDESTAL control (front panel) was intended to compensate only for the dc offset in the signal source, and for this reason, the range is only  $\pm 0.25V$ . Circuit modification is required for a larger range and is described in Sect. 8.7.

#### 6.1.1 Pedestal Adjustment with Signal-Source Offset Control

If the signal source has an offset control, begin the adjustment by replacing the signal source with a 50-ohm terminator (connected to the TC 307 INPUT). Connect a dc millivolt meter or a dc-coupled oscilloscope to the OUTPUT of the TC 307. **With** the gate signal removed (or below threshold), trim the PEDESTAL control so that no change in TC 307 output voltage occurs when the MODE switch is flipped. Note that the absolute voltage may be different from zero when this condition is attained. If this offset is excessive for the intended application, the internal offset controls must be readjusted. See (Sect. 6.2)

Remove the terminator and reconnect the signal source. Now trim **its** off set control as described above as if it

were the PEDESTAL adjustment. (If a meter is used for this adjustment, no pulses **must** emerge from the signal source because they will cause an error. The use of an oscilloscope permits pulses to be present because the baseline between pulses can be observed.)

### 6.1.2 Pedestal Adjustment, **Signal Source Offset not Adjustable**

Connect a dc millivolt meter or oscilloscope to the **OUTPUT** of the TC 307 and connect the signal source to the **INPUT**. Do not gate the TC 307 and do not allow pulses to issue from the signal source. Adjust the PEDESTAL control so that no change *in* dc output level occurs when the MODE switch is flipped.

If pulses cannot be turned off from the signal source, use the oscilloscope as the dc indicator to allow observation of the baseline between pulses. If the dc offset from the signal source is large enough to be outside the range of adjustment of the PEDESTAL control ( $\pm 0.25V$ ), the range *may* be increased by *an* internal modification. See Sect. 8.7.

### 6.2 Internal Adjustments

The internal adjustments are one-turn **screwdriver-**adjustable trimmers for normally closed zero (**N.C. ZERO**), normally open **zero** (**N.O. ZERO**), **OUTPUT ZERO\*** and thermal balance (**THERM BAL**). The trimmers are mounted on the circuit board. The left-side module cover plate must be removed for access. The adjustments can be made with the TC 307 mounted in a NIM bin if the TC 307 is inserted in the slot nearest the right side of the bin and if the modules to the left of the TC 307 are removed for access to the trimmers.

**THERM BAL** readjustment should never be required unless the input integrated circuit is replaced. Even then, readjustment may be unnecessary.

As in the **TBERM BAL** case, it is unlikely that the zero adjustments will require changes. However, an offset may result from power supply errors. The supply voltages should first be trimmed to within  $\pm 0.1V$  of their nominal values to confirm that **TC 307** readjustment is necessary.

---

\*Serial numbers above 209.



In addition to the TC 307 and its power source, a square-wave generator, an oscilloscope, and a 50 ohm terminator will be required.

Connect the system as shown in Fig. 6-2.




Fig. 6-2 Test setup for adjustment of TC 307 internal controls.

Set controls as follows:

OSCILLOSCOPE

Channel A	0.5V/Div dc coupled
Channel B	10 mV/Div dc coupled
sweep	1 ms/Div
Vert mode	CHOP

GENERATOR

Repetition rate	250 Hz
Amplitude	0.9V
Polarity	+

TC 307

MODE	N.C.
------	------

Trigger the oscilloscope sweep from Channel A (the generator channel) or directly from the generator. With **Channel B** disconnected, (or its INPUT switch set to GROUND), adjust the oscilloscope positioning controls to produce the waveforms shown in Fig. 6-3.




Fig. 6-3 Waveforms for TC 307 setup

Reconnect the TC 307 OUTPUT to Channel B. Adjust the N.O. ZERO and/or the OUTPUT ZERO (the two adjustments are virtually interchangeable) on the circuit **board** to zero output **volts**.

Switch the TC 307 MODE to N.O. and adjust the PEDESTAL control **for** zero volts output. If the control is out of range (which **may** occur if **IC1** or **IC2** are changed), alter the relative adjustments of the N.O. ZERO and OUTPUT ZERO and start over.

With the PEDESTAL and N.O. ZERO controls correctly set, the trace will not **move** from zero when the TC 307 MODE switch is flipped back and forth.

Increase the generator output to **1.5V**, causing the TC 307 gate to operate. Adjust the N.C. ZERO trimmer so that the Channel B trace remains at zero for both half-cycles **of** the square wave (the MODE switch setting is immaterial).

Next, adjust the THERM BAL trimmer so that the positive and negative transients **at the square-wave** transitions are symmetrical. See **Fig. 6-4**. Note that changing the setting of the **THERM** BAL trimmer will move the TC 307 **baseline**. After adjusting the trimmer for symmetrical transients, readjust the N.O. and N.C. **ZEROs**. This completes the alignment.




Fig. 6-4 THERM BAL trimmer adjustment

## 7.0 CIRCUIT OPERATION

A circuit diagram is at the end of this manual. The circuit is divided into four sections: Input Stage, Gate Switch, Output Stage, and Counting Monitor.

### 7.1. Input Stage

The input group is built around an **integrated** dual differential array (**IC1a** and **IC1b**). Depending on the state of the gate switch and the gate signal applied to pins 3 and 11 (Test Points 17 and 18) of **IC1**, **one** or the other half of the array will be conducting and will constitute the active input stage for the remainder of the inverting operational amplifier **Q1** through **Q5**. In the absence of gate signal, **IC1a** is normally conducting. The INPUT signal is routed through the MODE switch either to **IC1a** (N.O. position) or **IC1b** (N.C. position).

The input resistor is **R1** plus **R5** or **R6**, depending on the position of the MODE switch. **R8** and **R9** limit the loop gain of the input group for stability against oscillation.

The PEDESTAL bias is derived from **RO** and is injected through **R2** into the junction between **R1** and **R5**(**R6**). The diode network **D1-D4** limits the input signal to prevent any of the linear stages from saturating.

The **THERM BAL** control is used to set the operating voltage on the input **IC** transistors. Adjustable offset voltages are injected into the bases of these same transistors.

**R1, R5 (R6), and R10 (R11)** determine the gain of the input group. The gain is precisely 0.49 referred to a signal source with a 50 ohm internal resistance.

**Q1, Q2, and Q4** are gigahertz transistors.

## 7.2 Gate Switch

This stage is a Schmitt trigger. The loading on **Q6** and **Q7** was designed to give a symmetrical output such that the average voltage at pins 4 and 12 of **IC1** does not change when the gate is triggered. This symmetry minimizes gate transients. The switching voltage at pins 4 and 12 is set by **R58** and is just large enough to cause complete current switching between **IC1a** and **IC1b**. **D23-D26** limit the switching voltage of **Q6-Q7** and also speeds up the transition.

**D20-D22** protect **Q6-Q7** from excessive input voltage.

## 7.3 Output Stage

This is an inverting op amp built around **IC2**. **D9, D10,** and **D13** limit the base currents to **Q12** and **Q13** when the output is short circuited.

**IC3** and **IC4** are **15V** regulators which limit at 100 mA minimum; 200 mA typical. With the **OUTPUT IMPEDANCE** switch in the 50 ohm position, rated maximum output will be obtained if one of the two outputs is loaded with 50 ohms and may or may not be obtained if both outputs are loaded, depending on the saturation level of **IC3** or **IC4**. If the switch is in the **<1 OHM** position and either of the outputs is loaded with 50 ohms, a maximum signal level of 5V can be obtained, and the maximum may be 10V depending on the saturation level of **IC3** or **IC4**.

## 7.4 Counting Monitor

This circuit consists of a pair of comparators **IC5a** and **IC5b** connected as a univibrator. **IC5a** responds to positive polarity signals and **IC5b** to negative. **R92** and **R94** set the threshold of **IC5a**, and **D15** and **D16** set

the threshold of IC56. C52 and R95 determine the recovery time ( $\approx 1$  ms) after triggering. Q13 and Q14 are normally nonconducting. When IC5 is in the triggered state, Q13 and Q14 turn on, illuminating the LED D19. R103 and C56 constitute a saturating differentiator. For low count rates, the current into the LED is high so that the presence of a single pulse is visible. As the pulse rate increases, the average current drops off to maintain approximately constant LED illumination. C54 is a filter to prevent LED pulse current from appearing in the output signal;

## 8.0 OPERATING NOTES

### 6.1 Gate Control

The gate will remain in a switched state for the duration of a greater-than-1.3V gate signal. The duration and time of occurrence of the gate signal can be controlled by a Delay and Gate Generator, such as a TENNELEC TC 410A (0.1 to 110 us range for both delay and duration).

### a.2 Output Impedance

The output impedance may be set to 50 ohms  $\pm 1\%$  or to less than 1 ohm by a slide switch which is mounted on the circuit board. The switch is in the 50 ohm position when shipped. Unless there is a particular need for the  $\ll 1$  ohm output impedance, it is recommended that the switch be kept in the 50 ohm position.

### 8.3 $\ll 1$ ohm Operation

In the  $\ll 1$  ohm switch position, output loads on the TC 307 of 1k produce at most a 0.1% drop in signal height. However, except for very short connecting cables, cable termination probably will be required to avoid ringing. Critical cable lengths, unless terminated, may cause oscillation.

If the cable is terminated in 50 ohms, a **10V** output signal **will** require the output stage to feed **200 mA** into the termination. Because the voltage regulators which supply the output stage are guaranteed for a saturation level of only **100 mA**, not all TC 307s will furnish a **10V** signal into a **50 ohm** load, (However, the TC **307** will not be damaged by such operations).

a.4 50-ohm Operation

In the 50 ohm position,, a 50 ohm **±1%** resistor is interposed between the output stage and the front panel **OUTPUT** connector, and a separate **50 ohm ±1%** resistor is interposed between the output stage and the rear panel connector. The two output connectors are effectively isolated. from each other by the low common **resistance** of the output amplifier.

In the. 50 ohm position, cable termination should not be required for output cable lengths up to or **greater** than 100 feet, and then only if **signal risetime** is less than approximately **0.2 usec.** Termination will halve the signal height.

With an internal impedance of 50 ohms, 1k loading will cause signal attenuation **of** 5%. If two instruments with the same input resistance are driven **from a** TC 307, signal loss can be halved by driving one from the front-panel output connector and the other from the rear, but this arrangement usually **is unnecessary.**

If the output is loaded with 50 ohms (terminated cable in the **50-ohm** output impedance setting), the signal will be halved. Linearity will not be affected up to **5V** output in this instance **if** only one output connector is loaded, but if both connectors are loaded with **50 ohms**, the unit **may** not be able to **drive** the outputs to full voltage.

8.5 Placement in a System

The usual placement of a linear gate in a pulse height measurement system is between the output of the shaping amplifier and the following pulse height discriminator or multichannel analyzer. In situations where there is a high background of unwanted pulses and the desired Pulses are of relatively low rate and are time correlated with the gating signal, it is tempting to place the gate between the preamplifier and shaping

amplifier (keeping the gate in the normally closed mode) in an attempt to reduce pileup in the shaping amplifier. This is not a recommended arrangement for the following reasons:

1. Preamplifiers usually have a decay time which is many times that of the pulse width existing in the shaping amplifier. The result is that even at moderate count rates, the preamplifier output rarely **has an** opportunity to fully recover between pulses. **If the** gate is opened for a desired pulse, the **instantaneous** preamplifier output voltage will constitute a pedestal which is superimposed on the desired pulse. This pedestal will vary in amplitude and polarity from instant to instant, making for a very noisy measurement.

2. To avoid the preceding problem, the user may be tempted to reduce the decay time of the preamplifier, **in effect** placing the **system's** first differentiator **there**. This arrangement is undesirable for these reasons:

2a. Moving the first differentiator to the preamplifier can quadruple the system noise level.

2b. Superimposed on the preceding noise will be an added component which results from a gating operation in the part **of** the system which precedes the low pass portion of the shaping network.

2c. The equivalent **input** noise of the TC 307 is 3X to 7X greater than that **of** the usual shaping amplifier. This does not **affect** the system performance when the gate follows the main amplifier, but it may if the gate **precedes it**.

2d. The dc output stability of preamplifiers is poor. A changing dc level will result in a changing pedestal **in** the gate, further adding to the measured noise level.

## 6.6

### AC vs DC Coupling

A linear gate should never be ac coupled (capacitively coupled) to an amplifier because baseline shifts directly proportional to the count rate will occur. **Any** change in baseline will appear in the gate as a pedestal.

If the output of the gate is ac coupled to a following instrument, a baseline shift proportional to count rate also will occur unless the following instrument contains a baseline restorer.

0.7 Increasing the Range of the PEDESTAL Control

The range of the PEDESTAL control can be approximately doubled (from  $\pm 0.25V$  to  $\pm 0.5V$ ) by shunting R2 with a 20-kilohm resistor. For a further increase in range, the 50k PEDESTAL control should be replaced with a 10k, 2-watt unit and R2 further reduced. If R2 is reduced to 2 kilohms, the range will be increased to  $\pm 3V$ . However, this change will also reduce the gain of the TC 307 to a value which varies from 0.92 with the control at its midrange setting to 0.84 with the control at one of its end settings.

9.0 SHIPPING DAMAGE

Upon receipt of the instrument, examine it for shipping damage. Damage claims should be filed with the carrier. The claims agent should receive a full report: a copy of that report should be sent to TENNELEC, Inc., P.O. Box 2560, Oak Ridge, Tennessee 37830-2560. The model number and serial number of the instrument must be included in the report. Any remedial action taken by TENNELEC, Inc., will be based on the information contained in this report.

10.0 SERVICING

In the event of a component failure, replacement may be done in the field or the instrument may be returned to our plant for repair. There will be no charge for repairs that fall within the warranty.

11.0 WARRANTY

In connection with TENNELEC's warranty (inside front cover), TENNELEC suggests that if a fault develops, the customer should immediately notify the TENNELEC Customer Service Manager. He may be able to prescribe repairs and send replacement parts which will enable YOU to get the instrument operating sooner and at less expense than if you returned it.