

Target changing procedure

Preliminaries:

- Stop run
- Call operator to insert Faraday cup RT-ECR
- Push red button (user lockout beam blocker)
- Insert wall-plug and beam-stop for the S3 vault (on the control PC)
- Unplug Big Brother
- Turn off object scintillator voltage if needed
- Enable **REMOTELY** with the program I250StepMtr.chx
- I250StepMtr.chx : Enable 1250X-R (target; bad encoder), 1251Y-R (MCP0), 1250Y-R (MCP1)
- **On PC set the new target position (see note at PC) Program is called I250StepMtr.chx**

When finished with moving target:

- Disable **REMOTELY** with the program I250StepMtr.chx
- Don't forget to re-plug big brother

